

Anti-IL-10 antibody ab9969

★★★★★ [4 Abreviews](#) [77 References](#) [画像数 1](#)

製品の概要

製品名	Anti-IL-10 antibody
製品の詳細	Rabbit polyclonal to IL-10
由来種	Rabbit
アプリケーション	適用あり: ELISA, WB, ICC/IF, Sandwich ELISA
種交差性	交差種: Mouse, Rat
免疫原	Recombinant fragment corresponding to Rat IL-10 aa 19-178. (Peptide available as ab9970)
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Lyophilized:Reconstitute with 200µl of sterile water. Please note that if you receive this product in liquid form it has already been reconstituted as described and no further reconstitution is necessary.
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	No preservative, sterile filtered
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
軽鎖の種類	unknown

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab9969の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Can be blocked with <u>Recombinant rat IL-10 protein (Active) (ab9970)</u> .
ICC/IF	★★★★★ (1)	Use at an assay dependent concentration.
Sandwich ELISA		Use a concentration of 0.5 - 2 µg/ml.

ターゲット情報

機能	Inhibits the synthesis of a number of cytokines, including IFN-gamma, IL-2, IL-3, TNF and GM-CSF produced by activated macrophages and by helper T-cells.
組織特異性	Produced by a variety of cell lines, including T-cells, macrophages, mast cells and other cell types.
配列類似性	Belongs to the IL-10 family.
細胞内局在	Secreted.

画像



Sandwich ELISA ab9969. Allows the detection of at least 0.2 - 0.4 ng/well of recombinant Rat IL-10.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors