

Anti-IL-11 antibody [MM0361-9T27] ab89887

[1 References](#) [画像数 1](#)

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-IL-11 antibody [MM0361-9T27] |
| 製品の詳細 | Mouse monoclonal [MM0361-9T27] to IL-11 |
| 由来種 | Mouse |
| アプリケーション | 適用あり: WB, Neutralising |
| 種交差性 | 交差種: Human |
| 免疫原 | Recombinant full length protein corresponding to Human IL-11. |
| ポジティブ・コントロール | WB: Human placenta tissue lysate |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|----------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | Constituent: PBS |
| 精製度 | Protein G purified |
| 特記事項(精製) | The IgG fraction of culture supernatant was purified by Protein G affinity chromatography and 0.2 µm filtered. |
| ポリモノ | モノクローナル |
| クローン名 | MM0361-9T27 |
| アイソタイプ | IgG2 |

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab89887の使用に適用されず**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|--------------|-----------|---|
| WB | | 1/500 - 1/2000. Predicted molecular weight: 21 kDa. |
| Neutralising | | Use at an assay dependent concentration. |

ターゲット情報

| | |
|-------|---|
| 機能 | Directly stimulates the proliferation of hematopoietic stem cells and megakaryocyte progenitor cells and induces megakaryocyte maturation resulting in increased platelet production. |
| 配列類似性 | Belongs to the IL-6 superfamily. |
| 細胞内局在 | Secreted. |

画像



Anti-IL-11 antibody [MM0361-9T27] (ab89887) at 1/500 dilution +
Human placenta tissue lysate

Predicted band size: 21 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors