

### Anti-IL-12 p35 antibody ab66462

#### 製品の概要

---

|          |   |
|----------|---|
| 製品名      | Anti-IL-12 p35 antibody   |
| 製品の詳細    | Rabbit polyclonal to IL-12 p35  |
| 由来種      | Rabbit  |
| 特異性      | This antibody reacts specifically with the 22kDa IL-12 p35  |
| アプリケーション | <b>適用あり:</b> IHC-Fr, WB   |
| 種交差性     | <b>交差種:</b> Pufferfish  |
| 免疫原      | Synthetic peptide corresponding to IL-12 p35 (C terminal).  |
| 特記事項     | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

#### 製品の特性

---

|        |   |
|--------|---|
| 製品の状態  | Liquid  |
| 保存方法   | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー  | Constituent: Whole serum  |
| 精製度    | Whole antiserum   |
| ポリ/モノ  | ポリクローナル   |
| アイソタイプ | IgG   |

#### アプリケーション

---

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab66462の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項                                     |
|----------|-----------|--|
| IHC-Fr   |           | Use at an assay dependent concentration. |
| WB       |           | Use at an assay dependent concentration. |

## ターゲット情報

|              |   |
|--------------|---|
| <b>機能</b>    | Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated Killer cells, and stimulate the production of IFN-gamma by resting PBMC. |
| <b>配列類似性</b> | Belongs to the IL-6 superfamily.  |
| <b>細胞内局在</b> | Secreted.   |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors