

# Anti-IL-20R2 antibody [20RNTC] ab95824

**5 References** [画像数 2](#)

### 製品の概要

製品名	Anti-IL-20R2 antibody [20RNTC]
製品の詳細	Rat monoclonal [20RNTC] to IL-20R2
由来種	Rat
特異性	ab95824 recognizes Mouse and Human IL-20R2 in either of its two heterodimeric complexes formed by pairing with IL-20RA or IL-22RA1
アプリケーション	<b>適用あり:</b> Flow Cyt, ICC/IF
種交差性	<b>交差種:</b> Mouse, Human
免疫原	Recombinant fragment corresponding to Mouse IL-20R2 (extracellular).
ポジティブ・コントロール	Fixed and permeabilized A375 cells
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituent: 1% PBS</p>
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	20RNTC
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション


The Abpromise guarantee      Abpromise保証は、次のテスト済みアプリケーションにおけるab95824の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use a concentration of 1 - 10 µg/ml. <b>ab18450</b> - Rat monoclonal IgG2a, is suitable for use as an isotype control with this antibody.
ICC/IF		Use a concentration of 1 - 10 µg/ml.


ターゲット情報

機能	The IL20RA/IL20RB dimer is a receptor for IL19, IL20 and IL24. The IL22RA1/IL20RB dimer is a receptor for IL20 and IL24.
組織特異性	Widely expressed with highest levels in skin and testis. Highly expressed in psoriatic skin.
配列類似性	Belongs to the type II cytokine receptor family. Contains 2 fibronectin type-III domains.
細胞内局在	Membrane.

画像



Immunofluorescent microscopy of fixed and permeabilized A375 cells using 10 µg/ml of ab95824 (right) or Rat IgG2a - Isotype Control Purified (left) followed by Anti-Rat IgG Biotin and then Streptavidin FITC.



Mouse IL-20R2-GFP 293 transfectants stained with 0.25ug  
ab95824 + anti-rat IgG PE

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors