


### Anti-IL-25 antibody ab108530

**2 References**   [画像数 2](#)

#### 製品の概要

|              |   |
|--------------|---|
| 製品名          | Anti-IL-25 antibody   |
| 製品の詳細        | Rabbit polyclonal to IL-25  |
| 由来種          | Rabbit  |
| アプリケーション     | <b>適用あり:</b> WB, IHC-P  |
| 種交差性         | <b>交差種:</b> Human<br><b>交差が予測される動物種:</b> Rat   |
| 免疫原          | Recombinant full length protein corresponding to Human IL-25 aa 33-177.<br>Database link: <a href="#">BC069565</a>  |
| ポジティブ・コントロール | Human fetal lung and kidney lysates; Human fetal prostate tissue.   |
| 特記事項         | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

#### 製品の特性

|        |   |
|--------|---|
| 製品の状態  | Liquid  |
| 保存方法   | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.          |
| バッファー  | pH: 7.20<br>Preservative: 0.02% Sodium azide<br>Constituent: PBS<br><br>Buffer contains non-immune IgG. |
| 精製度    | Immunogen affinity purified   |
| ポリ/モノ  | ポリクローナル   |
| アイソタイプ | IgG   |

アプリケーション


The Abpromise guarantee      Abpromise保証は、 次のテスト済みアプリケーションにおけるab108530の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項  |
|----------|-----------|---|
| WB       |           | 1/500 - 1/1000. Predicted molecular weight: 20 kDa. |
| IHC-P    |           | 1/100 - 1/500.                                      |

ターゲット情報

|       |  |
|-------|--|
| 機能    | Induces activation of NF-kappa-B and stimulates production of the proinflammatory chemokine IL-8. Proinflammatory cytokine favoring Th2-type immune responses. |
| 組織特異性 | Expressed at low levels in several tissues, including brain, kidney, lung, prostate, testis, spinal cord, adrenal gland, and trachea.                          |
| 配列類似性 | Belongs to the IL-17 family.   |
| 細胞内局在 | Secreted.  |

画像




All lanes : Anti-IL-25 antibody (ab108530) at 1/500 dilution

Lane 1 : Human fetal lung lysate

Lane 2 : Human fetal kidney lysate

Predicted band size: 20 kDa



ab108530 at 1/100 dilution staining IL 17E in Human fetal prostate by Immunohistochemistry, Formalin-fixed, Paraffin-embedded tissue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors