

Anti-Integrin alpha 4+beta 7 antibody [RM0059-2G18] ab73261

画像数 1

製品の概要

製品名	Anti-Integrin alpha 4+beta 7 antibody [RM0059-2G18]
製品の詳細	Rat monoclonal [RM0059-2G18] to Integrin alpha 4+beta 7
由来種	Rat
アプリケーション	適用あり: Flow Cyt, IP, IHC-P
種交差性	交差種: Mouse
免疫原	Cell lysates of mouse lymphocytes
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituent: PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	RM0059-2G18
アイソタイプ	IgG2a

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab73261の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
IHC-P		1/100 - 1/200.

ターゲット情報

関連性	Integrin alpha 4 beta 7, also known as, the lymphocyte Peyer's patch adhesion molecule (LPAM1) is a member of the integrin family of cell surface receptors. It is expressed primarily on mucosal lymphocytes, but is also present on NK cells and eosinophils. The alpha 4 beta 7 heterodimer mediates the binding of lymphocytes to its ligand, mucosal vascular addressin (MAdCAM1) on the high endothelial venules, thereby directing the homing of lymphocytes into Peyer's patches and intestinal lamina propria.
細胞内局在	Membrane; Single pass type I membrane protein.

画像



ab73261 at 1/200 dilution, staining Integrin alpha 4+beta 7 in LPS exposed mouse spleen tissue section by Immunohistochemistry (Formalin/PFA fixed paraffin-embedded sections).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Integrin alpha 4+beta 7 antibody [RM0059-2G18] (ab73261)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors