

Anti-Integrin alpha V + beta 6 antibody [10D5] - BSA and Azide free ab77906

★★★★☆ [6 Abreviews](#) [16 References](#)

製品の概要

製品名	Anti-Integrin alpha V + beta 6 antibody [10D5] - BSA and Azide free
製品の詳細	Mouse monoclonal [10D5] to Integrin alpha V + beta 6 - BSA and Azide free
由来種	Mouse
特異性	ab77906 recognises alpha V beta 6 and does not recognize alpha V or other alpha V integrins.
アプリケーション	適用あり: Flow Cyt, IP, Blocking 適用なし: IHC-P or WB
種交差性	交差種: Mouse, Rat, Human
免疫原	Tissue, cells or virus corresponding to Mouse Integrin alpha V + beta 6. Wild type murine keratinocytes.
ポジティブ・コントロール	H441 cells
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituent: PBS
キャリア・フリー	はい
精製度	Protein A purified
特記事項(精製)	Sterile filtered
ポリ/モノ	モノクローナル

クローン名	10D5
アイソタイプ	IgG2a

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab77906の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration. Flow cytometry should be done with buffers containing 1 mM Mg, since the epitope recognized by this antibody is divalent cation dependent. ab170191 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.
IP		Use at an assay dependent concentration. because this antibody is an IgG2a, protein G or rabbit anti mouse IgG is helpful for capture.
Blocking		Use at an assay dependent concentration. Inhibition of alpha V beta 6 mediated cell adhesion and cell migration are achieved at a concentration of 100 µg/mL using beta 6 transfected cells tested against TGF beta LAP.

追加情報 Is unsuitable for IHC-P or WB.

ターゲット情報

機能	The alpha-V integrins are receptors for vitronectin, cytotactin, fibronectin, fibrinogen, laminin, matrix metalloproteinase-2, osteopontin, osteomodulin, prothrombin, thrombospondin and vWF. They recognize the sequence R-G-D in a wide array of ligands. In case of HIV-1 infection, the interaction with extracellular viral Tat protein seems to enhance angiogenesis in Kaposi's sarcoma lesions.
配列類似性	Belongs to the integrin alpha chain family. Contains 7 FG-GAP repeats.
細胞内局在	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors