

Anti-Interferon alpha 2 antibody [N39] ab24101

2 References

製品の概要

製品名	Anti-Interferon alpha 2 antibody [N39]
製品の詳細	Mouse monoclonal [N39] to Interferon alpha 2
由来種	Mouse
特異性	Specific for Interferon alpha 2 and does not crossreact with Interferon alpha 1.
アプリケーション	適用あり: ELISA, Blocking, WB 適用なし: IHC-P
種交差性	交差種: Human
免疫原	E.coli-derived recombinant full length Interferon alpha 2c (Human)
エピトープ	Carboxyl-terminal domain residues 112 to 148.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.3 Preservative: 0.05% Sodium azide Constituent: PBS
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	N39
ミエローマ	NS1
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab24101の使用に適用されず
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
Blocking		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

追加情報 Is unsuitable for IHC-P.

ターゲット情報

機能 Produced by macrophages, IFN-alpha have antiviral activities.

配列類似性 Belongs to the alpha/beta interferon family.

細胞内局在 Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors