


### Anti-ISOC2 antibody ab118201

面格粉

医薬用外毒物

#### 製品の概要

製品名	Anti-ISOC2 antibody
製品の詳細	Rabbit polyclonal to ISOC2
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to amino acids 26-203 of Human ISOC2 (BC017344).
ポジティブ・コントロール	HepG2 and Human fetal liver lysates; Human fetal liver tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Lyophilized:Reconstitute in 200µl sterile water
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituent: 1% BSA
精製度	Immunogen affinity purified
特記事項(精製)	ab118201 is purified by a protein affinity column.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション

### The Abpromise guarantee

**Abpromise保証は、**次のテスト済みアプリケーションにおけるab118201の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 22 kDa (predicted molecular weight: 22 kDa).
IHC-P		1/100 - 1/500.

## ターゲット情報


### 配列類似性

Belongs to the isochorismatase family.

### 細胞内局在

Mitochondrion. Cytoplasm. Nucleus. Localizes to the nucleus in the presence of CDKN2A.

## 画像




**All lanes :** Anti-ISOC2 antibody (ab118201) at 1/500 dilution

**Lane 1 :** PANC-1 cell

**Lane 2 :** Human fetal liver lysate

**Predicted band size:** 22 kDa

**Observed band size:** 22 kDa



ab118201, at 1/100 dilution, staining ISOC2 in formalin-fixed, paraffin-embedded Human fetal liver tissue showing nuclear staining, by Immunohistochemistry.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors