

# Anti-JunB antibody [204C4a] ab53543

**1 References** [画像数 1](#)

### 製品の概要

---

製品名	Anti-JunB antibody [204C4a]
製品の詳細	Mouse monoclonal [204C4a] to JunB
由来種	Mouse
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Recombinant fragment
免疫原	Recombinant fragment of JunB (Human)
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.4 Preservative: 0.05% Sodium azide Constituents: 1% BSA, PBS
精製度	Protein G purified
特記事項 (精製)	Purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG depleted (approximately 95%) fetal bovine serum
ポリ/モノ	モノクローナル
クローン名	204C4a
アイソタイプ	IgG1

### アプリケーション

---


**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab53543の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Detects a band of approximately 36 kDa (predicted molecular weight: 36 kDa).

#### ターゲット情報

機能	Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3'.
配列類似性	Belongs to the bZIP family. Jun subfamily. Contains 1 bZIP domain.
翻訳後修飾	Ubiquitinated by ITCH, leading to its degradation.
細胞内局在	Nucleus.

#### 画像



Anti-JunB antibody [204C4a] (ab53543) + Immunized recombinant protein

**Predicted band size:** 36 kDa

**Observed band size:** 36 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors