

Anti-Kallikrein 5 antibody ab28565

★★★★☆ [1 Abreviews](#) [5 References](#) [画像数 1](#)

製品の概要

製品名	Anti-Kallikrein 5 antibody
製品の詳細	Rabbit polyclonal to Kallikrein 5
由来種	Rabbit
特異性	This antibody reacts with Kallikrein 5, but does not react with the other Kallikreins (1-7 and 9-15).
アプリケーション	適用あり: WB
種交差性	交差種: Recombinant fragment
免疫原	Synthetic peptide based on the Kallikrein loop area of human Kallikrein 5. Immunogen の所有権に関して (Peptide available as ab41423 .)
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.05% Sodium azide Constituent: 50% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab28565の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/5000. Predicted molecular weight: 32 kDa. 1/1000 (when using colorimetric substrates such as BCIP/NBT)- 1/5000 (chemiluminescent substrates). Dilution optimised using Chromogenic detection.

ターゲット情報

機能	May be involved in desquamation.
組織特異性	Expressed in skin, breast, brain and testis. Expressed at the stratum granulosum of palmar skin.
配列類似性	Belongs to the peptidase S1 family. Kallikrein subfamily. Contains 1 peptidase S1 domain.
細胞内局在	Secreted.

画像



All lanes : Anti-Kallikrein 5 antibody (ab28565) at 1/1000 dilution

Lane 1 : Kallikrein 5 at 0.05 µg

Lane 2 : Kallikrein 5 at 0.01 µg

Lane 3 : Kallikrein 5 at 0.001 µg

Predicted band size: 32 kDa

Observed band size: 38 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors