


Anti-KAT13C / NCOA2 antibody ab10491

★★★★☆ **4 Abreviews** **13 References** 画像数 3

製品の概要

製品名	Anti-KAT13C / NCOA2 antibody
製品の詳細	Rabbit polyclonal to KAT13C / NCOA2
由来種	Rabbit
特異性	We have data to indicate that this antibody may not cross react with <i>Xenopus laevis</i> . However, this has not been conclusively tested and expression levels may vary in certain cell lines/tissues.
アプリケーション	適用あり: WB, IHC-P, IP
種交差性	交差種: Mouse, Human 交差が予測される動物種: Rat, Rabbit, Guinea pig, Cow, Dog, Pig, Chimpanzee, Ferret, Gorilla, Chinese hamster 
免疫原	Synthetic peptide corresponding to Human KAT13C/ NCOA2 aa 700-800. Database link: Q15596
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7 Preservative: 0.1% Sodium azide Constituents: 0.021% PBS, 1.764% Sodium citrate, 1.815% Tris
精製度	Immunogen affinity purified
特記事項(精製)	Antibodies were affinity purified using the peptide immobilized on solid support.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab10491の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★☆ (3)	1/1000 - 1/10000. Detects a band of approximately 170 kDa (predicted molecular weight: 175 kDa).
IHC-P		Use at an assay dependent concentration.
IP	★★★★☆ (1)	Use 6µg for 10 cells.

ターゲット情報

機能	Transcriptional coactivator for steroid receptors and nuclear receptors. Coactivator of the steroid binding domain (AF-2) but not of the modulating N-terminal domain (AF-1). Required with NCOA1 to control energy balance between white and brown adipose tissues.
関連疾患	Note=Chromosomal aberrations involving NCOA2 may be a cause of acute myeloid leukemias. Inversion inv(8)(p11;q13) generates the KAT6A-NCOA2 oncogene, which consists of the N-terminal part of KAT6A and the C-terminal part of NCOA2/TIF2. KAT6A-NCOA2 binds to CREBBP and disrupts its function in transcription activation.
配列類似性	Belongs to the SRC/p160 nuclear receptor coactivator family. Contains 1 bHLH (basic helix-loop-helix) domain. Contains 1 PAS (PER-ARNT-SIM) domain.
ドメイン	Contains four Leu-Xaa-Xaa-Leu-Leu (LXXLL) motifs. The LXXLL motifs are essential for the association with nuclear receptors and are, at least in part, functionally redundant. The LLXXLXXXL motif is involved in transcriptional coactivation and CREBBP/CBP binding. Contains 2 C-terminal transcription activation domains (AD1 and AD2) that can function independently.
翻訳後修飾	Phosphorylated upon DNA damage, probably by ATM or ATR.
細胞内局在	Nucleus.

画像



All lanes : Anti-KAT13C / NCOA2 antibody (ab10491) at 1/1000 dilution

Lane 1 : Wild-type HEK-293T cell lysate

Lane 2 : NCOA2 knockout HEK-293T cell lysate

Lane 3 : PC-3 cell lysate

Lane 4 : HeLa cell lysate


Lysates/proteins at 20 µg per lane.

Performed under reducing conditions.

Predicted band size: 175 kDa

Observed band size: 160 kDa

False colour image of Western blot: Anti-KAT13C / NCOA2 antibody staining at 1/1000 dilution, shown in green; Mouse anti-Alpha Tubulin [DM1A] ([ab7291](#)) loading control staining at 1/20000 dilution, shown in red. In Western blot, ab10491 was shown to bind specifically to KAT13C / NCOA2. A band was observed at 160 kDa in wild-type HEK-293T cell lysates with no signal observed at this size in NCOA2 knockout cell line [ab265068](#) (knockout cell lysate [ab258530](#)). To generate this image, wild-type and NCOA2 knockout HEK-293T cell lysates were analysed. First, samples were run on an SDS-PAGE gel then transferred onto a nitrocellulose membrane. Membranes were blocked in 5% milk in TBS-0.1 % Tween® 20 (TBS-T) before incubation with primary antibodies overnight at 4°C. Blots were washed four times in TBS-T, incubated with secondary antibodies for 1 h at room temperature, washed again four times then imaged. Secondary antibodies used were Goat anti-Rabbit IgG H&L (IRDye® 800CW) preabsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed ([ab216776](#)) at 1/20000 dilution.



Immunoprecipitation - Anti-KAT13C / NCOA2 antibody (ab10491)

ab10491 used for WB at 0.1 µg/ml (A)

Lane 1: Whole cell lysate HeLa, 50 ug

Lane 2: Whole cell lysate HeLa, 15 ug

Lane 3: Whole cell lysate HEK293T, 50 ug


Lane 4: Whole cell lysate Jurkat, 50 ug

ab10491 was used at 1 µg/ml for WB and used for IP at 6 µg/mg lysate (B)

Lane 1-3: ab10491]

Lane4: Control IgG

Detection: Chemiluminescence with exposure times of 3 minutes (A) and 10 seconds (B).



Western blot - Anti-KAT13C / NCOA2 antibody (ab10491)

All lanes : Anti-KAT13C / NCOA2 antibody (ab10491) at 1 µg/ml

Lane 1 : TCMK-1 cells

Lane 2 : 4T1 cells

Lane 3 : CT26.WT cells

Lysates/proteins at 50 µg per lane.

Predicted band size: 175 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors