

Anti-KIF3A antibody ab11259

★★★★★ [11 Abreviews](#) [41 References](#) [画像数 4](#)

製品の概要

製品名	Anti-KIF3A antibody
製品の詳細	Rabbit polyclonal to KIF3A
由来種	Rabbit
アプリケーション	適用あり: ICC/IF, WB, IHC-P
種交差性	交差種: Mouse, Rat, Cow, Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to amino acids 682-701 of Mouse KIF3A. This sequence is identical in human.
特記事項	<p>Storage in frost-free freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.40 Preservative: 0.097% Sodium azide Constituent: 0.0268% PBS
精製度	IgG fraction
特記事項(精製)	Whole antiserum is fractionated and then further purified by ion-exchange chromatography to provide the IgG fraction of antiserum that is essentially free of other rabbit serum proteins.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab11259の使用に適用されます**


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ICC/IF	★★★★★ (1)	Use at an assay dependent concentration.
WB	★★★★★ (9)	Use at an assay dependent concentration. Predicted molecular weight: 80 kDa. This concentration is determined by blotting using a rat brain homogenate and a whole cell extract of rat pheochromocytoma PC12 cell line. Staining of KIF3A is specifically inhibited by the KIF3A immunizing peptide.
IHC-P		1/100.

ターゲット情報


機能	Microtubule-based anterograde translocator for membranous organelles. Plus end-directed microtubule sliding activity in vitro.
配列類似性	Belongs to the kinesin-like protein family. Kinesin II subfamily. Contains 1 kinesin-motor domain.
細胞内局在	Cytoplasm > cytoskeleton.

画像




Paraffin-embedded rat brain tissue stained for KIF3A with ab11259 at 0.1 µg/mL in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-KIF3A antibody (ab11259)




Paraffin-embedded rat brain tissue stained for KIF3A with ab11259 at 0.1 $\mu\text{g}/\text{mL}$ in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-KIF3A antibody (ab11259)



Paraffin-embedded mouse brain tissue stained for KIF3A with ab11259 at 0.4 $\mu\text{g}/\text{mL}$ in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-KIF3A antibody (ab11259)



Paraffin-embedded mouse brain tissue stained for KIF3A with ab11259 at 0.4 µg/mL in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-KIF3A antibody (ab11259)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors