

Anti-KREMEN1 antibody [RM0120-6K45] ab86636

1 References [画像数 1](#)

製品の概要

製品名	Anti-KREMEN1 antibody [RM0120-6K45]
製品の詳細	Rat monoclonal [RM0120-6K45] to KREMEN1
由来種	Rat
特異性	This antibody detects mouse KREMEN1 in Western blotting.
アプリケーション	適用あり: WB
種交差性	交差種: Mouse
免疫原	Fusion protein corresponding to Mouse KREMEN1 aa 1-400.
特記事項	<p>This antibody was produced from a hybridoma (mouse myeloma fused with spleen cells from a rat immunized with mouse recombinant protein of the extracellular domain of KREMEN1).</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Lyophilized:Reconstitute the antibody with 200 ul sterile PBS for a final concentration of 0.5 mg/ml
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituent: PBS
精製度	Protein A purified
特記事項 (精製)	Purified from culture supernatant and lyophilized from a 0.2 µm filtered solution in PBS.
ポリ/モノ	モノクローナル
クローン名	RM0120-6K45
アイソタイプ	IgG2


アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab86636の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Predicted molecular weight: 52 kDa.

ターゲット情報

画像



Anti-KREMEN1 antibody [RM0120-6K45] (ab86636) + Lysates of mouse liver

Predicted band size: 52 kDa

Observed band size: 50 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors