

Anti-L Ficolin antibody [16] ab112558

製品の概要

製品名	Anti-L Ficolin antibody [16]
製品の詳細	Mouse monoclonal [16] to L Ficolin
由来種	Mouse
特異性	No cross reactivity is seen with M Ficolin (Ficolin 1) and H Ficolin (Ficolin 3).
アプリケーション	適用あり: Sandwich ELISA, WB
種交差性	交差種: Human
免疫原	Purified prokaryote expressed recombinant Human L Ficolin and full-length L Ficolin (Ficolin 2) protein adsorbed onto aluminum hydroxide gel.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.1% Sodium azide Constituents: 2.9% Sodium chloride, 97% Phosphate Buffer
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	16
ミエローマ	Sp2
アイソタイプ	IgG1
軽鎖の種類	kappa

