

### Anti-L Glutamate antibody ab9440

★★★★★ 5 Abreviews 13 References 画像数 1

**医薬用外劇物**

#### 製品の概要

<b>製品名</b>	Anti-L Glutamate antibody
<b>製品の詳細</b>	Rabbit polyclonal to L Glutamate
<b>由来種</b>	Rabbit
<b>特異性</b>	The antibody is calibrated against a spectrum of antigens to assure hapten selectivity. No measurable cross-reactivity (<1:1000) was detected against glutamate in peptides or proteins. Fixed tissue cross-reactivity was tested with known targets at the recommended dilution. No measurable glutaraldehyde-fixed tissue cross-reactivity (<1:1000) was detected against L-alanine, gamma-aminobutyrate, D/L-aspartate, agmatine, guanidine, D/L-arginine, L-citrulline, L-cysteine, D/L-glutamine, glutathione, glycine, L-lysine, L-ornithine, L-serine, taurine, L-threonine, L-tryptophan, L-tyrosine. Modest cross-reactivity (1:20) was detected against D-glutamate. Significant cross-reactivity (1:8) was detected against free NAAG in competition assays (NAAG is not retained in aldehyde-fixed tissue).
<b>アプリケーション</b>	<b>適用あり:</b> IHC-FoFr, IHC-Fr
<b>種交差性</b>	<b>交差種:</b> Species independent
<b>免疫原</b>	Chemical/ Small Molecule corresponding to L Glutamate (Glutaraldehyde).
<b>特記事項</b>	<p>The product is optimized for HPI/EHPI with gold or fluorescence detection using etched plastic sections. Filter diluted reagents with 0.2 mm syringe filters before use on EM grids. Enzyme-linked visualizations can be used but will compress the signal dynamic range and are less sensitive.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

<b>製品の状態</b>	Liquid
<b>保存方法</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

バッファー	Preservative: 0.05% Thimerosal (merthiolate) Constituents: Phosphate Buffer, 1% Goat Serum
精製度	IgG fraction
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab9440の使用に適用されます**


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-FoFr	★★★★☆ (1)	Use at an assay dependent concentration. Use with frozen or vibratome sections is possible but will not yield optimal images as IgGs penetrate aldehyde cross-linked tissue poorly and most amino acids are present at such high levels that prozone effects occur.
IHC-Fr	★★★★☆ (1)	Use at an assay dependent concentration.

## ターゲット情報

**関連性** L Glutamate is a non-essential amino acid naturally occurring in the L-form. It is the most common excitatory neurotransmitter in the Central Nervous System.

## 画像



Immunohistochemistry (Frozen sections) - Anti-L Glutamate antibody (ab9440)

This image is courtesy of an Abreview submitted by Patricia Bogdanov

ab9440 staining L Glutamate in Mouse retina tissue sections by Immunohistochemistry (IHC-Fr - frozen sections). Tissue was fixed with methanol, permeabilized with 0.025% TBS-Triton X-100 and blocked with 5% serum for 1 hour. Samples were incubated with primary antibody (1/100 in PBS + 1% BSA) for 18 hours at 4°C. An Alexa Fluor® 594-conjugated Goat anti-rabbit IgG polyclonal (1/200) was used as the secondary antibody.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors