

Anti-Lass5 antibody ab63927

画像 1

医薬用外劇物

製品の概要

| | |
|--------------|--|
| 製品名 | Anti-Lass5 antibody |
| 製品の詳細 | Rabbit polyclonal to Lass5 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Human |
| 免疫原 | Synthetic peptide surrounding amino acid 186 of human Lass5 |
| ポジティブ・コントロール | Jurkate cell lysate |
| 特記事項 | <p>This product is manufactured by BioVision, an Abcam company and was previously called 3830 LASS5 Antibody. 3830-100 is the same size as the 100 µg size of ab63927.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| バッファー | pH: 7.20 Preservative: 0.01% Thimerosal (merthiolate) Constituents: PBS, 30% Glycerol (glycerin, glycerine), 0.5% BSA |
| 精製度 | Immunogen affinity purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab63927の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| WB | | Use a concentration of 0.5 - 4 µg/ml. Detects a band of approximately 25 kDa (predicted molecular weight: 46 kDa). |

ターゲット情報

| | |
|-------|---|
| 機能 | May be either a bona fide (dihydro)ceramide synthase or a modulator of its activity. When overexpressed in cells is involved in the production of sphingolipids containing mainly one fatty acid donor (N-linked palmitoyl- (C16) ceramide) in a fumonisin B1-independent manner. |
| 配列類似性 | Contains 1 homeobox DNA-binding domain. Contains 1 TLC (TRAM/LAG1/CLN8) domain. |
| 細胞内局在 | Nucleus membrane. Endoplasmic reticulum membrane. |

画像



Anti-Lass5 antibody (ab63927) at 0.5 µg/ml + Jurkat cell lysate

Predicted band size: 46 kDa

Observed band size: 25 kDa

Additional bands at: 40 kDa (possible isoform), 50 kDa (possible isoform)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors