

Anti-LCAT antibody [EPR1384Y] ab51060

リコンビナント **RabMAb**

5 References **画像数 6**

製品の概要

製品名	Anti-LCAT antibody [EPR1384Y]
製品の詳細	Rabbit monoclonal [EPR1384Y] to LCAT
由来種	Rabbit
アプリケーション	適用あり: IP, WB, IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide within Human LCAT (C terminal). The exact sequence is proprietary. Database link: <u>P04180</u>
ポジティブ・コントロール	WB: Human serum and plasma lysates IHC-P: Human brain and liver tissue IP: Human plasma
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information <u>see here</u>.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <u>RabMAb[®] patents</u>.</p> <p>Mouse, Rat: We have preliminary internal testing data to indicate this antibody may not react with these species. Please contact us for more information.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.2</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 0.05% BSA, 40% Glycerol (glycerin, glycerine), 59% PBS</p>
精製度	Protein A purified
ポリ/モノ	モノクローナル

クローン名	EPR1384Y
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab51060の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご相談ください。

アプリケーション	Abreviews	特記事項
IP		1/20. For unpurified use at 1/100.
WB		1/10000. Detects a band of approximately 62 kDa (predicted molecular weight: 50 kDa).
IHC-P		<p>1/100 - 1/250. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. See <u>IHC antigen retrieval protocols.</u></p> <p>The immunostaining was performed on a Leica Biosystems BOND® RX instrument.</p>

ターゲット情報

機能 Central enzyme in the extracellular metabolism of plasma lipoproteins. Synthesized mainly in the liver and secreted into plasma where it converts cholesterol and phosphatidylcholines (lecithins) to cholesteryl esters and lysophosphatidylcholines on the surface of high and low density lipoproteins (HDLs and LDLs). The cholesterol ester is then transported back to the liver. Has a preference for plasma 16:0-18:2 or 18:0-18:2 phosphatidylcholines. Also produced in the brain by primary astrocytes, and esterifies free cholesterol on nascent APOE-containing lipoproteins secreted from glia and influences cerebral spinal fluid (CSF) APOE- and APOA1 levels. Together with APOE and the cholesterol transporter ABCA1, plays a key role in the maturation of glial-derived, nascent lipoproteins. Required for remodeling high-density lipoprotein particles into their spherical forms.

組織特異性	Expressed mainly in brain, liver and testes. Secreted into plasma and cerebral spinal fluid. In liver, expressed in HEPG2 hepatocytes.
--------------	--

関連疾患

Defects in LCAT are the cause of lecithin-cholesterol acyltransferase deficiency (LCATD) [MIM:245900]; also called Norum disease. LCATD is a disorder of lipoprotein metabolism characterized by inadequate esterification of plasmatic cholesterol. Two clinical forms are recognized: familial LCAT deficiency and fish-eye disease. Familial LCAT deficiency is associated with a complete absence of alpha and beta LCAT activities and results in esterification anomalies involving both HDL (alpha-LCAT activity) and LDL (beta-LCAT activity). It causes a typical triad of diffuse corneal opacities, target cell hemolytic anemia, and proteinuria with renal failure.

Defects in LCAT are a cause of fish-eye disease (FED) [MIM:136120]; also known as dyslipoproteinemic corneal dystrophy or alpha-LCAT deficiency. FED is due to a partial LCAT deficiency that affects only alpha-LCAT activity. It is characterized by low plasma HDL and corneal opacities due to accumulation of cholesterol deposits in the cornea ('fish-eye').

配列類似性 Belongs to the AB hydrolase superfamily. Lipase family.


翻訳後修飾

O- and N-glycosylated. O-glycosylation on Thr-431 and Ser-433 consists of sialylated galactose beta 1-->3N-acetylgalactosamine structures. N-glycosylated sites contain sialylated triantennary and/or biantennary complex structures.

細胞内局在

Secreted. Secreted into blood plasma. Produced in astrocytes and secreted into cerebral spinal fluid.

画像



Anti-LCAT antibody [EPR1384Y] (ab51060) at 1/10000 dilution (Purified) + Human plasma lysates at 15 µg


Secondary

Goat Anti-Rabbit IgG (HRP) with minimal cross-reactivity with human IgG at 1/2000 dilution

Predicted band size: 50 kDa


Observed band size: 62 kDa

Blocking/Diluting Buffer and concentration: 5% NFDM/TBST



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human liver tissue sections labeling LCAT with purified ab51060 at 1/100 dilution (1.03 µg/mL). Heat mediated antigen retrieval was performed using Bond™ Epitope Retrieval Solution 2 (pH 9.0) . Rabbit specific IHC polymer detection kit HRP/DAB ([ab209101](#)) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain.

The immunostaining was performed on a Leica Biosystems BOND® RX instrument.



Western blot - Anti-LCAT antibody [EPR1384Y]
(ab51060)


Anti-LCAT antibody [EPR1384Y] (ab51060) at 1/10000 dilution
(unpurified) + Human serum at 10 µg

Secondary

Goat anti-Rabbit HRP at 1/2000 dilution

Predicted band size: 50 kDa

Observed band size: 62 kDa



Immunoprecipitation - Anti-LCAT antibody
[EPR1384Y] (ab51060)

ab51060 (Purified) at 1/20 dilution (0.5ug) immunoprecipitating
LCAT in Human plasma.

Lane 1 (input): Human plasma 10ug

Lane 2 (+): ab51060 & Human plasma


Lane 3 (-): Rabbit monoclonal IgG (**ab172730**) instead of ab51060
in Human plasma

For western blotting, VeriBlot for IP Detection Reagent (HRP)

(**ab131366**) was used at 1/1000 dilution.

Blocking and diluting buffer: 5% NFDm/TBST.

50 kDa band in input lane should be human IgG heavy chain.



Immunohistochemistry (Formalin/PFA-fixed paraffin-
embedded sections) - Anti-LCAT antibody
[EPR1384Y] (ab51060)

ab51060 (unpurified) at 1/100 dilution staining LCAT in paraffin
embedded human brain tissue by Immunohistochemistry. Heat
mediated antigen retrieval was performed before commencing with
IHC staining protocol.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-LCAT antibody [EPR1384Y] (ab51060)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors