

Anti-LTBR antibody [5G11] ab65089

8 References [画像数 1](#)

製品の概要

製品名	Anti-LTBR antibody [5G11]
製品の詳細	Rat monoclonal [5G11] to LTBR
由来種	Rat
アプリケーション	適用あり: Flow Cyt 適用なし: IHC-Fr or IHC-P
種交差性	交差種: Mouse
免疫原	Fusion protein corresponding to LTBR. Please refer to the reference below.
特記事項	<p>This antibody is an agonist of LTBR in vitro and is provided in a low endotoxin format.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at +4°C.
バッファー	Constituents: PBS, 0.1% BSA
精製度	Protein G purified
特記事項(精製)	ab65089 was 0.2µm filtered.
一次抗体 備考	This antibody is an agonist of LTBR in vivo.
ポリ/モノ	モノクローナル
クローン名	5G11
アイソタイプ	IgG2a

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab65089の使用に適用されず**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		1/50. ab18450 - Rat monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

追加情報 Is unsuitable for IHC-Fr or IHC-P.


ターゲット情報

機能 Receptor for the heterotrimeric lymphotoxin containing LTA and LTB, and for TNFS14/LIGHT. Promotes apoptosis via TRAF3 and TRAF5. May play a role in the development of lymphoid organs.

配列類似性 Contains 4 TNFR-Cys repeats.

細胞内局在 Membrane.

画像



Detection of LT beta receptor in Wehi-167 cells. Red, black and blue line represent the isotype control, cells only and ab65089 at 10 µg/ml, respectively.

Flow Cytometry - Anti-LTBR antibody [5G11]
 (ab65089)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors