


### Anti-Lubricin/MSF antibody ab28484

★★★★★ [7 Abreviews](#) [50 References](#) [画像数 1](#)

#### 製品の概要

<b>製品名</b>	Anti-Lubricin/MSF antibody
<b>製品の詳細</b>	Rabbit polyclonal to Lubricin/MSF
<b>由来種</b>	Rabbit
<b>特異性</b>	ab28484 recognises Lubricin/MSF.
<b>アプリケーション</b>	<b>適用あり:</b> IHC-P, ICC/IF, WB
<b>種交差性</b>	<b>交差種:</b> Mouse, Human, African green monkey <b>交差が予測される動物種:</b> Horse 
<b>免疫原</b>	Synthetic peptide corresponding to Human Lubricin/MSF aa 1356-1374. Sequence: CLPNIRKQPDGYDYAFSKDQ <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
<b>特記事項</b>	The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.  If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 製品の特性

<b>製品の状態</b>	Liquid
<b>保存方法</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>バッファー</b>	Preservative: 0.05% Sodium azide Constituents: PBS, 0.1% BSA
<b>精製度</b>	Immunogen affinity purified
<b>ポリモノ</b>	ポリクローナル
<b>アイソタイプ</b>	IgG

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab28484の使用に適用されず**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。


アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (4)	Use at an assay dependent concentration.
ICC/IF	★★★★☆ (1)	Use at an assay dependent concentration.
WB	★★★★☆ (1)	Use a concentration of 2 µg/ml. Detects a band of approximately 280 kDa (predicted molecular weight: 156 kDa).

### ターゲット情報

**関連性** 28484 is a glycoprotein encoded by the PRG4 gene. It has been identified as a component of joint lubrication and synovial homeostasis by reducing friction in joints and preventing synovial cell overgrowth. It is a large, water soluble molecule that acts as a carrier for insoluble surface-active phospholipid (SAPL). Depletion of 28484 function has been associated with camptodactyly-arthropathy-coxa vara-pericarditis syndrome (CACP), an arthritic-like autosomal recessive disorder.

**細胞内局在** Secreted

### 画像



Anti-Lubricin/MSF antibody (ab28484) at 2 µg/ml + non-reduced human synovial fluid at 3 µg

**Predicted band size:** 156 kDa

Primary antibody was diluted in antibody diluent using a Fast Western Blot kit and followed by the optimized HRP reagent. Chemiluminescent detection was performed using SuperSignal West Dura.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors