

### Anti-LYPLAL1 antibody ab77640

[2 References](#) [画像数 1](#)

#### 製品の概要

製品名	Anti-LYPLAL1 antibody
製品の詳細	Goat polyclonal to LYPLAL1
由来種	Goat
アプリケーション	<b>適用あり:</b> WB
種交差性	交差種: Human 交差が予測される動物種: Rat, Horse, Cow, Dog, Chimpanzee, Macaque monkey, Orangutan

免疫原	Synthetic peptide: DRFKITNDCPEHLE , corresponding to internal sequence amino acids 79-92 of Human LYPLAL1 (NP_620149.1) <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
-----	--

**特記事項**

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Store in the dark. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 0.5% BSA, Tris buffered saline
精製度	Immunogen affinity purified
特記事項 (精製)	ab77640 was purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリ/モノ	ポリクローナル

## アプリケーション

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab77640の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		

## 追加情報

Peptide ELISA: Use at an assay dependent dilution. Antibody detection limit dilution: 1/16000.  
WB: Use at a concentration of 0.3 - 1 µg/ml. Detects a band of approximately 26 kDa (predicted molecular weight: 26 kDa).

Not yet tested in other applications.  
Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

**配列類似性** Belongs to the AB hydrolase 2 family.

**細胞内局在** Cytoplasm.

## 画像



Anti-LYPLAL1 antibody (ab77640) at 0.3 µg/ml + Human adipose lysate in RIPA buffer at 35 µg

**Predicted band size:** 26 kDa

**Observed band size:** 28 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors