


Anti-MAGP2 antibody ab47826

[1 References](#) [画像数 1](#)

製品の概要

製品名	Anti-MAGP2 antibody
製品の詳細	Rabbit polyclonal to MAGP2
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Chimpanzee 
免疫原	Synthetic peptide corresponding to a region near the carboxy terminal end of MAGP2 protein (Human).
ポジティブ・コントロール	MAGP-2 transfected HEK293 lysate
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab47826の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		

追加情報


ELISA: 1/5000 - 1/20000.
 IF (microscopy): Use at an assay dependent dilution.
 WB: 1/250 - 1/1000. Predicted molecular weight: 19 kDa. .

 Not yet tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能 Component of the elastin-associated microfibrils.
配列類似性 Belongs to the MFAP family.
翻訳後修飾 Forms intermolecular disulfide bonds either with other MAGP-2 molecules or with other components of the microfibrils.
細胞内局在 Secreted > extracellular space > extracellular matrix.

画像



All lanes : Anti-MAGP2 antibody (ab47826) at 1/100 dilution

Lane 1 : Supernatants from non-transfected cells

Lane 2 : Secreted MAGP-2

Lane 3 : MAGP-2 transfected HEK293 lysate

Predicted band size: 19 kDa

Observed band size: 19 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors