


# Anti-MAP2 antibody - Neuronal Marker ab32454

★★★★★ [35 Abreviews](#) [256 References](#) [画像数 7](#)

### 製品の概要

製品名	Anti-MAP2 antibody - Neuronal Marker
製品の詳細	Rabbit polyclonal to MAP2 - Neuronal Marker
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P, WB, ICC/IF
種交差性	<b>交差種:</b> Mouse, Rat, Human <b>交差が予測される動物種:</b> Goat, Cat, Lizard 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p> <p>Batches of this product that have a concentration &lt; 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help.</p>
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション


**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab32454の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (8)	Use a concentration of 1 µg/ml. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.
WB	★★★★★ (1)	Use at an assay dependent concentration. Predicted molecular weight: 199 kDa.
ICC/IF	★★★★★ (11)	1/1000.

## ターゲット情報

機能	The exact function of MAP2 is unknown but MAPs may stabilize the microtubules against depolymerization. They also seem to have a stiffening effect on microtubules.
配列類似性	Contains 3 Tau/MAP repeats.
翻訳後修飾	Phosphorylated at serine residues in K-X-G-S motifs by MAP/microtubule affinity-regulating kinase (MARK1 or MARK2), causing detachment from microtubules, and their disassembly (By similarity). Isoform 2 is probably phosphorylated by PKA at Ser-323, Ser-354 and Ser-386 and by FYN at Tyr-67.
細胞内局在	Cytoplasm, cytoskeleton.

## 画像




ab32454 staining MAP2 in Rat Primary Neurons DIV14 cells. The cells were fixed with 100% methanol (5 min), permeabilized with 0.1% PBS-Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1%PBS-Tween for 1h. The cells were then incubated overnight at 4°C with ab32454 at 1µg/ml and **ab7291**, Mouse monoclonal [DM1A] to alpha Tubulin - Loading Control. Cells were then incubated with **ab150081**, Goat polyclonal Secondary Antibody to Rabbit IgG - H&L (Alexa Fluor® 488), pre-adsorbed at 1/1000 dilution (shown in green) and **ab150120**, Goat polyclonal Secondary Antibody to Mouse IgG - H&L (Alexa Fluor® 594), pre-adsorbed at 1/1000 dilution (shown in pseudocolour magenta). Nuclear DNA was labelled with DAPI (shown in blue).

Also suitable in cells fixed with 4% paraformaldehyde (10 min).

Image was acquired with a high-content analyser (Operetta CLS, Perkin Elmer) and a maximum intensity projection of confocal


sections is shown.



Immunocytochemistry/ Immunofluorescence - Anti-MAP2 antibody - Neuronal Marker (ab32454)

Müller A et al. Monitoring Astrocytic Proteome Dynamics by Cell Type-Specific Protein Labeling. PLoS One 10:e0145451 (2015).

Rat astrocytes (DIV 22) cells were fixed with 4% paraformaldehyde and stained for MAP2 (Green) using ab32454 at 1/1000 dilution in ICC/IF analysis.



Western blot - Anti-MAP2 antibody - Neuronal Marker (ab32454)

**All lanes** : Anti-MAP2 antibody - Neuronal Marker (ab32454) at 1 µg/ml

**Lane 1** : Mouse Brain (day 0) Tissue Lysate

**Lane 2** : Mouse Brain Tissue Lysate

**Lane 3** : Rat Brain Tissue Lysate

Lysates/proteins at 20 µg per lane.

### Secondary

**All lanes** : Goat polyclonal to Rabbit IgG - H&L - Pre-Adsorbed (HRP) at 1/50000 dilution


**Predicted band size:** 199 kDa

**Observed band size:** 268,280 kDa

**Exposure time:** 2 minutes

Blocking buffer: 2% BSA


Gel type: TA



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MAP2 antibody - Neuronal Marker (ab32454)

This image is courtesy of an Abreview by Carl Hobbs.

ab32454 staining MAP2 in mouse brain tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with formaldehyde and blocked with 1% BSA for 10 minutes at 21°C; antigen retrieval was by heat mediation in a citrate buffer. Samples were incubated with primary antibody (1/3000 in blocking buffer) for 2 hours at 21°C in TBS/BSA/azide. An undiluted Biotin-conjugated goat anti-rabbit IgG polyclonal was used as the secondary antibody.



Western blot - Anti-MAP2 antibody - Neuronal Marker (ab32454)


Gonzales PK et al. Transcriptome analysis of genetically matched human induced pluripotent stem cells disomic or trisomic for chromosome 21. PLoS One 13:e0194581 (2018).

**All lanes** : Anti-MAP2 antibody - Neuronal Marker (ab32454) at 1/1000 dilution

**Lanes 1-2** : Trisomic induced Pluripotent Stem Cells (iPSCs)

**Lanes 3-4** : Disomic induced Pluripotent Stem Cells (iPSCs)

**Predicted band size:** 199 kDa




Immunocytochemistry/ Immunofluorescence - Anti-MAP2 antibody - Neuronal Marker (ab32454)

This image is courtesy of an Abreview submitted by Babben Tinner

ab32454 staining MAP2 in rat cryopreserved embryonic cortical neurons cells by ICC/IF

(Immunocytochemistry/immunofluorescence). Cells were fixed with paraformaldehyde with picric acid. Samples were incubated with primary antibody (1/2000 in 10mM PBS + 0.3% Triton-X) for 12 hours at 4°C. An Alexa Fluor® 488-conjugated donkey anti-rabbit IgG polyclonal (1/200) was used as the secondary antibody.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MAP2 antibody - Neuronal Marker (ab32454)

IHC image of MAP2 staining in human cerebral cortex FFPE section, performed on a Leica Bond™ system using the standard protocol F. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH6, epitope retrieval solution 1) for 20 mins. The section was then incubated with ab32454, 1 µg/ml, for 15 mins at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors