

### Anti-MCM5 antibody [CRCT5.1 (A2.7A3)] ab6154

★★★★★ [1 Abreviews](#) [6 References](#) [画像数 1](#)

#### 製品の概要

製品名	Anti-MCM5 antibody [CRCT5.1 (A2.7A3)]
製品の詳細	Mouse monoclonal [CRCT5.1 (A2.7A3)] to MCM5
由来種	Mouse
特異性	Recognises human Mcm-5
アプリケーション	<b>適用あり:</b> IHC-P
種交差性	<b>交差種:</b> Human
免疫原	Recombinant full length protein corresponding to Human MCM5.
ポジティブ・コントロール	Cervical Carcinoma Tissue
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium azide
精製度	Protein G purified
特記事項 (精製)	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	CRCT5.1 (A2.7A3)
ミエローマ	Sp2
アイソタイプ	IgG2b

## アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab6154の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/25 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## ターゲット情報

配列類似性	Belongs to the MCM family. Contains 1 MCM domain.
細胞内局在	Nucleus.

## 画像



ab6154 staining MCM5 in human tonsil by Immunohistochemistry  
(Formalin/PFA-fixed paraffin-embedded sections).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MCM5 antibody  
[CRCT5.1 (A2.7A3)] (ab6154)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors