

Anti-MD1 antibody [MD113] ab11860

製品の概要

製品名	Anti-MD1 antibody [MD113]
製品の詳細	Rat monoclonal [MD113] to MD1
由来種	Rat
アプリケーション	適用あり: Flow Cyt 適用なし: IP or WB
種交差性	交差種: Mouse
免疫原	Full length native protein (purified). This information is proprietary to Abcam and/or its suppliers.
特記事項	<p>Prior to use in biological assays, antibody should be filter sterilised and dialysed to remove azide.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium azide Constituents: PBS, 0.1% BSA
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	MD113
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab11860の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration. ab37371 - Rat polyclonal IgG, is suitable for use as an isotype control with this antibody.

追加情報 Is unsuitable for IP or WB.

ターゲット情報

機能	May cooperate with CD180 and TLR4 to mediate the innate immune response to bacterial lipopolysaccharide (LPS) and cytokine production. Important for efficient CD180 cell surface expression.
組織特異性	Highly expressed in B-cells, monocytes and tonsil.
細胞内局在	Secreted > extracellular space. Associated with CD180 at the cell surface.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors