

Anti-Metapneumovirus Fusion Protein antibody [hMPV 24] ab94800

3 References

製品の概要

製品名	Anti-Metapneumovirus Fusion Protein antibody [hMPV 24]
製品の詳細	Mouse monoclonal [hMPV 24] to Metapneumovirus Fusion Protein
由来種	Mouse
特異性	ab94800 is useful for the identification and location of expression of the fusion glycoprotein of Human metapneumovirus (hMPV) of both types A and B.
アプリケーション	適用あり: ICC/IF, ELISA
種交差性	交差種: Metapneumovirus
免疫原	Tissue, cells or virus corresponding to Metapneumovirus Fusion Protein. hMPV strain NCL-145
特記事項	<p>Fusion partner: PS-NS/1-Ag4 Recommended growth conditions: Growth medium RPMI with 25mM Hepes, 2mM additional glutamine and 10% heat inactivated foetal calf serum</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	hMPV 24
アイソタイプ	IgG2a

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab94800の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 20 µg/ml.
ELISA		Use at an assay dependent dilution.

ターゲット情報

関連性 Fusion glycoprotein precursor belonging to the paramyxoviruses fusion glycoprotein family

細胞内局在 Membrane; Single-pass type I membrane protein

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors