

# Anti-Metapneumovirus Phosphoprotein antibody [hMPV 33] ab94801

## 3 References

### 製品の概要

|              |  |
|--------------|--|
| 製品名          | Anti-Metapneumovirus Phosphoprotein antibody [hMPV 33]   |
| 製品の詳細        | Mouse monoclonal [hMPV 33] to Metapneumovirus Phosphoprotein   |
| 由来種          | Mouse  |
| 特異性          | Identification and location of expression of human MPV phosphoprotein of both types A and B.   |
| アプリケーション     | <b>適用あり:</b> ELISA, ICC/IF   |
| 種交差性         | <b>交差種:</b> Metapneumovirus  |
| 免疫原          | Tissue, cells or virus corresponding to Metapneumovirus Phosphoprotein. Human hMPV strain NCL-174  |
| ポジティブ・コントロール | Human MPV infected LLC-MK2 cells   |
| 特記事項         | <p>Fusion partner: PS-NS/1-Ag4 Growth medium RPMI with 25mM Hepes, 2mM additional glutamine and 10% heat inactivated foetal calf serum</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

### 製品の特性

|       |   |
|-------|---|
| 製品の状態 | Liquid  |
| 保存方法  | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | Preservative: 0.02% Sodium azide<br>Constituent: 99.98% PBS   |
| 精製度   | Protein A purified  |
| ポリ/モノ | モノクローナル   |
| クローン名 | hMPV 33   |

ミエローマ P3-NS1/1-Ag4-1  
アイソタイプ IgG2a

## アプリケーション

**The Abpromise guarantee** Abpromise保証は、 次のテスト済みアプリケーションにおけるab94801の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項   |
|----------|-----------|--|
| ELISA    |           | Use at an assay dependent concentration.               |
| ICC/IF   |           | Use at an assay dependent concentration. Acetone fixed |

## ターゲット情報

**関連性** Human metapneumovirus (hMPV) is a member of the Paramyxoviridae family belonging to the Metapneumovirus genus.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors