

Anti-Mumps nucleoprotein antibody [7B10] ab9880

12 References [画像数 1](#)

製品の概要

製品名	Anti-Mumps nucleoprotein antibody [7B10]
製品の詳細	Mouse monoclonal [7B10] to Mumps nucleoprotein
由来種	Mouse
アプリケーション	適用あり: WB, Indirect ELISA
種交差性	交差種: Mumps virus
免疫原	Tissue, cells or virus corresponding to Mumps nucleoprotein. Recombinant mumps NP. Gloucester strain (Glouc 1/UK96)
特記事項	<p>This product was changed from ascites to tissue culture supernatant on 28/11/2017. Lot numbers higher than GR293636-4 and GR137414-5 will be from tissue culture supernatant. Please note that the dilutions may need to be adjusted accordingly.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.1% Sodium azide
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	7B10
ミエローマ	NS0/1
アイソタイプ	IgG2a

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab9880の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/5000. Predicted molecular weight: 66 kDa.
Indirect ELISA		1/1000 - 1/10000.

ターゲット情報

画像



Immunoblot of yeast-derived mumps NP using monoclonal antibody 7B10.

Lane 1 - prestained standard proteins

Lane 2 - mumps NP

Antibody 7B10 at a concentration of 0.5 µg/ml.

Immunoblot of yeast-derived mumps NP using monoclonal antibody 7B10.

Lane 1 - prestained standard proteins

Lane 2 - mumps NP

Antibody 7B10 at a concentration of 0.5 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors