

Anti-MYEOV antibody ab121387

[1 References](#) [画像数 2](#)

製品の概要

製品名	Anti-MYEOV antibody
製品の詳細	Rabbit polyclonal to MYEOV
由来種	Rabbit
アプリケーション	適用あり: ICC/IF, IHC-P
種交差性	交差種: Human
免疫原	antigen, corresponding to amino acids 10-149 of Human MYEOV.
ポジティブ・コントロール	IHC: Human cerebral cortex tissue. ICC/IF: Human cell line A-431
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab121387の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報

発生段階


Overexpressed in tumor cells lines with a t(11;14)(q13;q32) translocation.

画像



Immunocytochemistry/immunofluorescence analysis of human cell line A-431 cells labelling MYEOV with ab121387 at 2 µg/mL.

Immunocytochemistry/ Immunofluorescence - Anti-MYEOV antibody (ab121387)



ab121387 at 1/400 dilution staining MYEOV in Paraffin-embedded Human cerebral cortex tissue by Immunohistochemistry, showing strong cytoplasmic positivity in neuronal cells.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MYEOV antibody (ab121387)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors