


Anti-NAP1L1 antibody ab21630

★★★★★ [4 Abreviews](#) [5 References](#) [画像数 4](#)

製品の概要

製品名	Anti-NAP1L1 antibody
製品の詳細	Rabbit polyclonal to NAP1L1
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB, ICC/IF
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Dog, Xenopus laevis 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	Recombinant Human NAP1L1 protein (ab117213) can be used as a positive control in WB. This antibody gave a positive signal in the following whole cell lysates: U2OS (+/- ionizing radiation); HeLa; A431; HEK293. This antibody gave a positive signal in the following tissues: Formalin Fixed Paraffin Embedded Normal Human Skin.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p> <p>Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help.</p>

精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab21630の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 5 µg/ml. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.
WB	★★★★★ (3)	1/250 - 1/1000. Detects a band of approximately 52 kDa (predicted molecular weight: 45 kDa).
ICC/IF	★★★★★ (1)	Use a concentration of 1 µg/ml.

ターゲット情報

機能	May be involved in modulating chromatin formation and contribute to regulation of cell proliferation.
組織特異性	Ubiquitously expressed.
配列類似性	Belongs to the nucleosome assembly protein (NAP) family.
ドメイン	The acidic domains are probably involved in the interaction with histones.
翻訳後修飾	Polyglutamylated by TTL4, a modification that occurs exclusively on glutamate residues and results in polyglutamate chains on the gamma-carboxyl group. Some residues may also be monoglycylated but not polyglycylated due to the absence of functional TTL10 in human.
細胞内局在	Nucleus. Melanosome. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

画像



Lanes 1-2 : Anti-NAP1L1 antibody (ab21630) at 1/1000 dilution

Lanes 3-4 : Anti-NAP1L1 antibody (ab21630) at 1/500 dilution

Lanes 5-6 : Anti-NAP1L1 antibody (ab21630) at 1/250 dilution

Lanes 1 & 3 & 5 : U2OS cell lysate

Lanes 2 & 4 & 6 : U2OS cell lysate treated with ionizing radiation


Performed under reducing conditions.

Predicted band size: 45 kDa

Observed band size: 52 kDa

Additional bands at: 49 kDa (possible cleavage fragment), 49 kDa (possible cross reactivity)

For each lane, U2OS cells, either treated or not treated with ionizing radiation, were scraped from a 60 mm dish and added to 75µl of 2X Laemmli buffer. 20µl of these samples were loaded into each lane. ab21630 recognizes a major band of approximately 52 kDa corresponding closely in size to NAP1L1.



All lanes : Anti-NAP1L1 antibody (ab21630) at 1 µg/ml

Lane 1 : HeLa whole cell lysate at 20 µg

Lane 2 : A431 whole cell lysate at 20 µg

Lane 3 : HEK293 whole cell lysate at 20 µg

Lane 4 : HeLa whole cell lysate at 20 µg with Human NAP1L1 peptide (**ab22417**) at 1 µg/ml

Lane 5 : A431 whole cell lysate with Human NAP1L1 peptide (**ab22417**) at 1 µg/ml

Lane 6 : HEK293 whole cell lysate at 20 µg with Human NAP1L1 peptide (**ab22417**) at 1 µg/ml

Secondary


All lanes : Alexa fluor goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 45 kDa

Observed band size: 52 kDa

Additional bands at: 35 kDa (possible cross reactivity, but this band is not blocked), 49 kDa (possible cleavage fragment), 49 kDa (possible cross reactivity)


ab21630 recognizes a major band of approximately 52 kDa corresponding closely in size to NAP1L1. This band is competed away by the addition of the immunizing peptide, suggesting that this is a specific interaction.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NAP1L1 antibody (ab21630)

IHC image of ab21630 staining in human skin formalin fixed paraffin embedded tissue section, performed on a Leica Bond™ system using the standard protocol F. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH6, epitope retrieval solution 1) for 20 mins. The section was then incubated with ab21630, 5µg/ml, for 15 mins at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

For other IHC staining systems (automated and non-automated) customers should optimize variable parameters such as antigen retrieval conditions, primary antibody concentration and antibody incubation times.



Immunocytochemistry/ Immunofluorescence - Anti-NAP1L1 antibody (ab21630)

ICC/IF image of ab21630 stained A431 cells. The cells were 4% formaldehyde fixed (10 min) and then incubated in 1% BSA / 10% normal goat serum / 0.3M glycine in 0.1% PBS-Tween for 1h to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab21630, 5µg/ml) overnight at +4°C. The secondary antibody (green) was [ab96899](#), DyLight® 488 goat anti-rabbit IgG (H+L) used at a 1/250 dilution for 1h. Alexa Fluor® 594 WGA was used to label plasma membranes (red) at a 1/200 dilution for 1h. DAPI was used to stain the cell nuclei (blue) at a concentration of 1.43µM.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors