

Anti-Neuritin antibody ab64186

★★★★★ [1 Abreviews](#) [4 References](#) [画像数 2](#)

製品の概要

製品名	Anti-Neuritin antibody
製品の詳細	Rabbit polyclonal to Neuritin
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to Human Neuritin (internal sequence). 20 amino acid peptide from near the center of human Neuritin.
ポジティブ・コントロール	Human brain cortex tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.4</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p>
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab64186の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	Use at an assay dependent concentration. Predicted molecular weight: 15 kDa.
IHC-P		Use a concentration of 10 µg/ml.

ターゲット情報


機能	Promotes neurite outgrowth and especially branching of neuritic processes in primary hippocampal and cortical cells.
配列類似性	Belongs to the neuritin family.
細胞内局在	Cell membrane. Cell junction > synapse.

画像



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Neuritin antibody (ab64186)

ab64186 (10ug/ml) staining Neuritin in human brain, cortex tissue by immunohistochemistry using formalin-fixed, paraffin-embedded tissue.



Western blot - Anti-Neuritin antibody (ab64186)

Western blot analysis of neuritin in Daudi cell lysate using ab64186 at (A) 5 and (B) 10 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors