


# Anti-Neuronatin antibody ab27266

★★★★★ [3 Abreviews](#) [19 References](#) [画像数 3](#)

### 製品の概要

製品名	Anti-Neuronatin antibody
製品の詳細	Rabbit polyclonal to Neuronatin
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Mouse <b>交差が予測される動物種:</b> Rat, Human 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: 1% BSA, PBS</p> <p>Batches of this product that have a concentration &lt; 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help.</p>
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab27266の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 9 kDa (predicted molecular weight: 9 kDa). Abcam recommends using 1-3% Milk as the blocking agent.

## ターゲット情報

機能	May participate in the maintenance of segment identity in the hindbrain and pituitary development, and maturation or maintenance of the overall structure of the nervous system. May function as a regulatory subunit of ion channels.
配列類似性	Belongs to the neuronatin family.
発生段階	Abundant in 18-24 week old fetal brain. Postnatally its expression decline and only minimal levels were present in adulthood.

## 画像




Anti-Neuronatin antibody (ab27266) at 1 µg/ml + Mouse brain tissue lysate - total protein (0 days) (**ab7188**) at 20 µg

### Secondary

IRDye 680 Conjugated Goat Anti-Rabbit IgG (H+L) at 1/15000 dilution

Performed under reducing conditions.

**Predicted band size:** 9 kDa



Western blot - Anti-Neuronatin antibody (ab27266)

Anti-Neuronatin antibody (ab27266) at 1 µg/ml + Brain (Mouse)  
Tissue Lysate - normal tissue, 0 days old (blocked with 3% Milk) at  
50 µg

### Secondary

Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/10000 dilution


Developed using the ECL technique.

Performed under reducing conditions.

**Predicted band size:** 9 kDa

**Observed band size:** 9 kDa

**Exposure time:** 12 minutes



Western blot - Anti-Neuronatin antibody (ab27266)

Anti-Neuronatin antibody (ab27266) at 1 µg/ml + Mouse brain  
tissue lysate - total protein (0 days) ([ab7188](#)) at 10 µg

### Secondary

Goat polyclonal to Rabbit IgG - H&L - Pre-Adsorbed (HRP) at  
1/3000 dilution

Performed under reducing conditions.

**Predicted band size:** 9 kDa

**Observed band size:** 10 kDa

**Additional bands at:** 70-75 kDa. We are unsure as to the identity  
of these extra bands.

**Exposure time:** 5 minutes

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors