


Anti-Neuropeptide S antibody ab18252

★★★★☆ [2 Abreviews](#) [8 References](#) [画像数 2](#)

製品の概要

製品名	Anti-Neuropeptide S antibody
製品の詳細	Rabbit polyclonal to Neuropeptide S
由来種	Rabbit
特異性	This antibody is generated against the complete 20aa sequence of the fully processed rat NPS peptide, it should run at about 2kDa on WB. ab18252 detects mouse NPS peptide in IHC (PFA fixed) and in RIA in hypothalamic mouse brain tissue.
アプリケーション	適用あり: IHC (PFA fixed), RIA
種交差性	交差種: Rat 交差が予測される動物種: Mouse 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help. Immunogen affinity purified

ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab18252の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC (PFA fixed)		Use a concentration of 0.3 µg/ml.
RIA		Use a concentration of 0.1 - 0.3 µg/ml.

ターゲット情報

機能	Modulates arousal and anxiety. May play an important anorexigenic role (By similarity). Binds to its receptor NPSR1 with nanomolar affinity to increase intracellular calcium concentrations (PubMed:15312648, PubMed:16790440).
細胞内局在	Secreted.


画像



Immunohistochemistry (PFA fixed) - Anti-Neuropeptide S antibody (ab18252)

This image is courtesy of Sophie Pezet, Univ London Kings Coll, United Kingdom.

Immunofluorescent staining using the Neuropeptide S antibody (ab18252) in the rat hypothalamus. Picture taken using objectives X20 (A) or X40 (B). The staining is observed in thin and varicose fibers. This staining is blocked by preincubation with immunising peptides ([ab18698](#) and [ab18699](#)). ab18252 was used at 1/1000 (0.3µg/ml) incubated overnight at room temperature. Secondary antibody used was anti-rabbit Alexa Fluor 488 at 1/1000 incubated for 2 hours at room temperature. Rat brain tissue was perfusion fixed with 4% PFA followed by overnight post-fixation in the same fixative, cryoprotected in 20% sucrose and frozen in OCT. 30µm coronal sections were cut on a cyrostat and immunohistochemistry performed by the 'free floating' technique.



IC50: 46 fmoles; ab18252 used at a dilution of 1/3000; Tracer: [125I] Y10-hNPS, 12.000 cpm; Detection method: anti-rabbit SPA (Amersham).

Radioimmunoassay was performed in equilibrium, using the SPA technology (Amersham). 50µl antirabbit YSi SPA beads (RPN 140) were mixed with 50µl antiserum (1:3000), 50µl tracer (12,000 cpm) and with or without cold NPS peptide in increasing concentrations. All reagents were resuspended in 0.1M phosphate buffer, pH 7, 0.1% BSA. Incubation was for 48h at RT. Bound radioligand was detected by scintillation proximity in a beta counter (MicroBeta, Wallac).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors