


Anti-NFAT5 antibody ab110995

2 References [画像数 1](#)

製品の概要

製品名	Anti-NFAT5 antibody
製品の詳細	Rabbit polyclonal to NFAT5
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse 
免疫原	Synthesized peptide derived from the internal region of Human NFAT5.
ポジティブ・コントロール	Human breast carcinoma tissue
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), PBS, 0.88% Sodium chloride
精製度	Immunogen affinity purified
特記事項 (精製)	ab110995 is affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
ポリ/モノ	ポリクローナル
アイソタイプ	IgG


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab110995の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報

機能	Plays a role in the inducible expression of genes. Regulates hypertonicity-induced cellular accumulation of osmolytes.
組織特異性	Highest levels in skeletal muscle, brain, heart and peripheral blood leukocytes. Also expressed in placenta, lung, liver, kidney, pancreas, spleen, thymus, prostate, testis, ovary, small intestine and colon.
配列類似性	Contains 1 RHD (Rel-like) domain.
細胞内局在	Nucleus.

画像



ab110995, at a 1/50 dilution, staining NFAT5 in paraffin-embedded Human breast carcinoma tissue, in the absence (left) and presence (right) of the immunizing peptide.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors