

Anti-NLRC5 antibody ab105411

15 References [画像数 1](#)

製品の概要

製品名	Anti-NLRC5 antibody
製品の詳細	Rabbit polyclonal to NLRC5
由来種	Rabbit
特異性	ab105411 is predicted to not cross react with other NLRP family members.
アプリケーション	適用あり: WB
種交差性	交差種: Mouse
免疫原	An 18 amino acid peptide near the amino terminus of Human NLRC5 (NP_115582).
ポジティブ・コントロール	WB: EL4 cell lysate.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、 次のテスト済みアプリケーションにおけるab105411の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 205 kDa. Can be blocked with NLRC5 peptide (ab196047) .

ターゲット情報

機能	Probable regulator of the NF-kappa-B and type I interferon signaling pathways. May also regulate the type II interferon signaling pathway. Plays a role in homeostatic control of innate immunity and in antiviral defense mechanisms.
組織特異性	Expressed in spleen, thymus, lung, brain, tonsil, heart and prostate.
配列類似性	Belongs to the NLRP family. Contains 26 LRR (leucine-rich) repeats. Contains 1 NACHT domain.
細胞内局在	Cytoplasm.

画像



Anti-NLRC5 antibody (ab105411) at 1 µg/ml + EL4 cell lysate

Predicted band size: 205 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors