

Anti-NLRP13 antibody ab105410

画像数 3

製品の概要

製品名	Anti-NLRP13 antibody
製品の詳細	Rabbit polyclonal to NLRP13
由来種	Rabbit
特異性	ab105410 will recognize both isoforms of NLRP13. It is predicted to not crossreact with other NLRP family members.
アプリケーション	適用あり: ICC/IF, WB, ICC
種交差性	交差種: Human
免疫原	A 19 amino acid peptide near the N terminus of Human NLRP13.
ポジティブ・コントロール	K562 cell lysate
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	<p>pH: 7.2</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p>
精製度	Immunogen affinity purified
特記事項(精製)	Affinity chromatography purified via peptide column.
ポリモノ	ポリクローナル
アイソタイプ	IgG


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab105410の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 119 kDa.
ICC		Use at an assay dependent concentration.


ターゲット情報

機能	Involved in inflammation.
配列類似性	Belongs to the NLRP family. Contains 1 DAPIN domain. Contains 7 LRR (leucine-rich) repeats. Contains 1 NACHT domain.


画像



Immunofluorescence of NLRP13 in K562 cells with ab105410 at 20 ug/mL.



Immunocytochemistry of K562 cells with ab105410 at 2.5 µg/mL.



Anti-NLRP13 antibody (ab105410) at 1 µg/ml + K562 cell lysate

Predicted band size: 119 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors