

# Anti-NOTUM antibody ab106448

★★★★★ [1 Abreviews](#) [1 References](#) [画像数 2](#)

### 製品の概要

---

製品名	Anti-NOTUM antibody
製品の詳細	Rabbit polyclonal to NOTUM
由来種	Rabbit
アプリケーション	<b>適用あり:</b> ICC/IF, WB, IHC-P
種交差性	<b>交差種:</b> Mouse, Rat, Human
免疫原	A 19 amino acid peptide near the amino terminus of Human NOTUM (NP_848588).
ポジティブ・コントロール	WB: Human kidney tissue lysate IHC-P: Human brain tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

---


**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab106448の使用に適用されず  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.
WB	★☆☆☆☆ (1)	Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 56 kDa.
IHC-P		Use at an assay dependent concentration.

## ターゲット情報

機能	May deacetylate GlcNAc residues on cell surface glycans.
配列類似性	Belongs to the pectinacetyltransferase family.
細胞内局在	Secreted.

## 画像



Western blot - Anti-NOTUM antibody (ab106448)


**Lane 1** : Anti-NOTUM antibody (ab106448) at 0.5 µg/ml

**Lane 2** : Anti-NOTUM antibody (ab106448) at 1 µg/ml

**All lanes** : Human kidney tissue lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 56 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NOTUM antibody (ab106448)

ab106448 at 2.5µg/ml staining NOTUM in Human brain by Immunohistochemistry.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors