

Anti-NSP1 antibody [32D6] ab4641

★★★★★ [1 Abreviews](#) [8 References](#) [画像数 1](#)

製品の概要

製品名	Anti-NSP1 antibody [32D6]
製品の詳細	Mouse monoclonal [32D6] to NSP1
由来種	Mouse
特異性	This antibody was initially identified as directed against a nuclear pore complex antigen by immunofluorescence localization. A screen of a lgt11 expression library yielded a single positive clone carrying an insert bearing ~66% of the C-terminal portion of NSP1. To confirm that 32D6 possessed high affinity for NSP1, strain RS453, which expresses a shortened isoform of NSP1, was compared to the wildtype strain BJ5465. NSP1 in RS453 contains an internal deletion that removes the coding sequence for 6 FXFG repeats, which are not essential for function and encodes a protein that has been observed to migrate at approximately 85 kDa on SDS gels. The predicted size of wildtype NSP1 is 86.5 kDa, but NSP1 has been observed to migrate on SDS gels at ~100 kDa. In our gels, the apparent molecular mass of wildtype NSP1 was 108 kDa, whereas the short isoform of NSP1 in RS453 migrated at 91 kDa. The detection of two protein bands of apparent sizes 108 kDa and 91 kDa in BJ5465 and RS453, respectively, demonstrated that 32D6 recognized the pore complex protein NSP1.
アプリケーション	適用あり: ICC/IF, WB
種交差性	交差種: <i>Saccharomyces cerevisiae</i>
免疫原	Tissue, cells or virus corresponding to NSP1. Yeast nuclear preparations
特記事項	This antibody is a useful marker of yeast nuclear pores. The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

バッファー	Preservative: 0.065% Sodium azide Constituent: Tissue culture supernatant
精製度	Tissue culture supernatant
特記事項 (精製)	Sterile filtered.
一次抗体 備考	This antibody is a useful marker of yeast nuclear pores.
ポリ/モノ	モノクローナル
クローン名	32D6
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab4641の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF	★★★★★ (1)	1/1000. Yeast cells.
WB		1/10000 - 1/50000. Detects a band of approximately 108,91 kDa (predicted molecular weight: 85 kDa). 1/10000 (cell lysates) - 1/50000 (nuclear fractions). For other (non-ECL) detection methods, 1/1000 - 1/5000.

ターゲット情報

関連性	NSP1 functions as a component of the nuclear pore complex (NPC). It plays an important role in several nuclear transport pathways including poly(A)+ RNA, tRNA, pre-ribosome, signal recognition particle (SRP), and protein transport.
細胞内局在	Nuclear membrane: part of the nuclear pore complex, it presents a symmetric distribution with a cytoplasmic and a nucleoplasmic side.

画像



Immunocytochemistry/ Immunofluorescence - Anti-NSP1 antibody [32D6] (ab4641)

Image courtesy of an anonymous Abreview.

ab4641 (1/5000) detecting NSP1 in *Saccharomyces cerevisiae* DDY1810 (green). Cells were fixed in paraformaldehyde, permeabilized in methanol/acetone and counterstained with DAPI in order to highlight the nucleus (blue). For further experimental details please refer to Abreview.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors