

Anti-O-Linked N-Acetylglucosamine antibody [RL2] ab2739

★★★★★ **17 Abreviews** **162 References** 画像数 4

製品の概要

製品名	Anti-O-Linked N-Acetylglucosamine antibody [RL2]
製品の詳細	Mouse monoclonal [RL2] to O-Linked N-Acetylglucosamine
由来種	Mouse
アプリケーション	適用あり: ICC/IF, WB
種交差性	交差種: Rat, Human
免疫原	Tissue, cells or virus corresponding to O-Linked N-Acetylglucosamine. Specifically, isolated rat liver nuclear envelopes, which contain 8 O-Linked glycoproteins in the nuclear pore complex
ポジティブ・コントロール	ICC-IF: MCF7 cells. WB: Jurkat cells treated with 50 uM PugNAc; SH-SY5Y) whole cell lysate - treated with 50µM z-Pugnac; Rat Liver Nuclear Envelope lysate.
特記事項	<p>This antibody clone [RL2] is manufactured by Abcam.</p> <p>If you require this antibody in a particular buffer formulation or a particular conjugate for your experiments, please contact orders@abcam.com or you can find further information here.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: PBS, 6.97% L-Arginine</p>
精製度	Protein G purified
ポリ/モノ	モノクローナル

クローン名	RL2
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab2739の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF	★★★★★ (2)	Use a concentration of 5 - 10 µg/ml.
WB	★★★★★ (15)	Use a concentration of 1 µg/ml.

ターゲット情報


関連性 Many cellular proteins, including nuclear pore, oncogene, cytoskeletal, heat shock, viral and transcription regulatory proteins contain single O-linked N-acetylglucosamine (O-GlcNAc) residues attached to serine or threonine residues. It has been observed that O-GlcNAc glycosylated proteins tend to be under phosphorylated relative to unglycosylated proteins and that O-GlcNAc bearing proteins tend to be found in multimeric complexes. This has led to the suggestion that O-GlcNAc glycosylation may obscure phosphorylation sites and acts as a signaling mechanism or mediator of signaling.

画像



Immunocytochemistry/ Immunofluorescence - Anti-O-Linked N-Acetylglucosamine antibody [RL2] (ab2739)

ab2739 stained in MCF7 cells. Cells were fixed with 4% paraformaldehyde (10min) at room temperature and incubated with PBS containing 10% goat serum, 0.3 M glycine, 1% BSA and 0.1% triton for 1h at room temperature to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody ab2739 at 5µg/ml and [ab6046](#) (Rabbit polyclonal to beta Tubulin - Loading Control) at 1/1000 dilution overnight at +4°C. The secondary antibodies were [ab150080](#) (pseudo-colored red) and [ab150117](#) (colored green) used at 1 ug/ml for 1hour at room temperature. DAPI was used to stain the cell nuclei (colored blue) at a concentration of 1.43µM for 1hour at room temperature.



Western blot - Anti-O-Linked N-Acetylglucosamine antibody [RL2] (ab2739)

All lanes : Anti-O-Linked N-Acetylglucosamine antibody [RL2] (ab2739) at 1 µg/ml

Lanes 1 & 3 : Jurkat cells treated with 0 uM PugNAc

Lane 2 : Jurkat cells treated with 50 uM PugNAc (3 hours)

Lane 4 : Jurkat cells treated with 4 mM glucosamine and 50 uM PugNAc (3 hours)

Lysates/proteins at 20 µg per lane.

Secondary


All lanes : Goat Anti-Mouse IgG H&L (HRP) preadsorbed (**ab97040**) at 1/50000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 12 minutes

Jurkat cells were treated with either 50 uM PugNAc (**ab144670**) or 4 mM glucosamine + 50 uM PugNAc (**ab144670**) for three hours prior to harvest to stimulate O-linked glycosylation. The expected increase in glycosylation is observed in the treated lanes 2 & 4.



Western blot - Anti-O-Linked N-Acetylglucosamine antibody [RL2] (ab2739)

This image is courtesy of an anonymous Abreview

All lanes : Anti-O-Linked N-Acetylglucosamine antibody [RL2] (ab2739) at 1/3000 dilution

Lane 1 : Human neuroblastoma (SH-SY5Y) whole cell lysate - treated with 50µM z-Pugnac for 24 hours

Lane 2 : Human neuroblastoma (SH-SY5Y) whole cell lysate - untreated

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : HRP-conjugated horse anti-mouse IgG polyclonal

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 30 seconds



Anti-O-Linked N-Acetylglucosamine antibody [RL2] (ab2739) at 1 µg/ml + Rat Liver Nuclear Envelope at 10 µg

Secondary

Goat polyclonal to Mouse IgG - H&L - Pre-Adsorbed (HRP) at 1/3000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 1 minute

The antibody was tested against the immunogen (isolated rat liver nuclear envelopes, which contain 8 O-linked glycoproteins in the nuclear pore complex).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors