


Anti-Oligodendrocyte myelin glycoprotein antibody ab101567

面格粉 1

医薬用外毒物

製品の概要

製品名	Anti-Oligodendrocyte myelin glycoprotein antibody
製品の詳細	Rabbit polyclonal to Oligodendrocyte myelin glycoprotein
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to internal sequence amino acids 94-376 of Human Oligodendrocyte myelin glycoprotein (NM_002544).
ポジティブ・コントロール	Human fetal brain
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Lyophilized:Reconstitute in 200 ul sterile water.
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium azide
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab101567の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/100 - 1/500.

ターゲット情報

機能	Cell adhesion molecule contributing to the interactive process required for myelination in the central nervous system.
組織特異性	Oligodendrocytes and myelin of the central nervous system.
配列類似性	Contains 8 LRR (leucine-rich) repeats. Contains 1 LRRNT domain.
翻訳後修飾	O-glycosylated in its Ser/Thr-rich repeat domain.
細胞内局在	Cell membrane.

画像



ab101567, at a 1/100 dilution, showing membrane staining of Oligodendrocyte myelin glycoprotein in a formalin fixed, paraffin embedded section of Human fetal brain by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Oligodendrocyte myelin glycoprotein antibody (ab101567)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors