

Anti-OVGP1/OGP antibody ab74544

6 References [画像数 3](#)

製品の概要

製品名	Anti-OVGP1/OGP antibody
製品の詳細	Rabbit polyclonal to OVGP1/OGP
由来種	Rabbit
アプリケーション	適用あり: WB, ICC/IF
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to Human OVGP1/OGP (N terminal). Database link: NP_002548
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab74544の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 150 kDa.
ICC/IF		Use a concentration of 20 µg/ml.

ターゲット情報


機能	Binds to oocyte zona pellucida in vivo. May play a role in the fertilization process and/or early embryonic development.
組織特異性	Oviduct.
配列類似性	Belongs to the glycosyl hydrolase 18 family.
細胞内局在	Cytoplasmic vesicle > secretory vesicle. Secretory granules.

画像



PFA-fixed, Triton X-100 permeabilized HEK-293 (human epithelial cell line from embryonic kidney) cells stained for OVGP1/OGP (green) using ab74544 at 5 µg/ml in ICC/IF.

Immunocytochemistry/ Immunofluorescence - Anti-OVGP1/OGP antibody (ab74544)



Anti-OVGP1/OGP antibody (ab74544) at 1 µg/ml + human placenta tissue lysate

Observed band size: 150 kDa

Western blot - Anti-OVGP1/OGP antibody (ab74544)

Immunocytochemistry/ Immunofluorescence - Anti-OVGP1/OGP antibody (ab74544)

PFA-fixed, Triton X-100 permeabilized HEK-293 (human epithelial cell line from embryonic kidney) cells stained for OVGP1/OGP using ab74544 at 2.5 µg/ml in ICC/IF.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors