

# Anti-p15 INK4b antibody [DCS-114] ab44891

[1 Abreviews](#) [4 References](#)

## 製品の概要

製品名	Anti-p15 INK4b antibody [DCS-114]
製品の詳細	Mouse monoclonal [DCS-114] to p15 INK4b
由来種	Mouse
アプリケーション	<b>適用あり:</b> ELISA, ICC, IP, ICC/IF, IHC <b>適用なし:</b> WB
種交差性	<b>交差種:</b> Human <b>非交差種:</b> Mouse, Rat
免疫原	Recombinant full length protein corresponding to Human p15 INK4b.
ポジティブ・コントロール	Ls174T cells
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

## 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.1% Sodium azide Constituent: PBS
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	DCS-114
アイソタイプ	IgG1

## アプリケーション

**The Abpromise guarantee**      Abpromise保証は、 次のテスト済みアプリケーションにおけるab44891の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
ICC		Use at an assay dependent concentration.
IF		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration.
IHC		Use at an assay dependent concentration.

**追加情報**      Is unsuitable for WB.

## ターゲット情報

**機能**      Interacts strongly with CDK4 and CDK6. Potent inhibitor. Potential effector of TGF-beta induced cell cycle arrest.

**配列類似性**      Belongs to the CDKN2 cyclin-dependent kinase inhibitor family.  
Contains 4 ANK repeats.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors