


Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] ab32142

KO 評価済 リコンビナント RabMAb

[42 References](#) [画像数 9](#)

製品の概要

製品名	Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122]
製品の詳細	Rabbit monoclonal [Y122] to p38 beta/MAPK11 + p38 alpha/MAPK14
由来種	Rabbit
特異性	This antibody recognises p38. It is predicted to react with splice isoform CSBP1 according to sequence homology. Mouse cross reactivity has been tested by WB and IHC, Rat cross reactivity by WB only.
アプリケーション	適用あり: ICC/IF, WB, IHC-P 適用なし: Flow Cyt or IP
種交差性	交差種: Mouse, Human 交差が予測される動物種: Rat 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: HeLa, Jurkat, HEK-293T, K562 and MCF7 cell lysates. IHC-P: Human skin carcinoma. ICC/IF: MCF7 cells.
特記事項	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide

Constituents: 49% PBS, 50% Glycerol (glycerin, glycerine), 0.05% BSA

ポリモノ
クローン名
アイソタイプ

モノクローナル
Y122
IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab32142の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
ICC/IF		1/100 - 1/250.
WB		1/1000 - 1/10000. Detects a band of approximately 42 kDa (predicted molecular weight: 41 kDa).
IHC-P		Use at an assay dependent concentration. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

追加情報 Is unsuitable for Flow Cyt or IP.

ターゲット情報

細胞内局在 p38 alpha/MAPK14: Cytoplasm. Nucleus.

画像



All lanes : Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142) at 1/1000 dilution

Lane 1 : MAPK11 recombinant (**ab117219**)

Lane 2 : MAPK12 recombinant (**ab117221**)

Lane 3 : MAPK13 recombinant (**ab113869**)

Lane 4 : MAPK14 (p38) recombinant (**ab82188**)

Lysates/proteins at 0.5 µg per lane.


Performed under reducing conditions.

Predicted band size: 41 kDa

Observed band size: 43 kDa

Lanes 1 - 4: Green - ab32142 observed at 43 kDa.

ab32142 was shown to react with Anti-p38 antibody [Y122] in Western blot. Membranes were blocked in 100% Licor before incubation with ab32142 and overnight at 4 °C at a 1 in 1000 dilution. Blots were incubated with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preabsorbed (**ab216773**) secondary antibody at 1 in 20000 dilution for 1 h at room temperature before imaging.



Western blot - Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142)

All lanes : Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142) at 1/1000 dilution


Lane 1 : Recombinant Human p38 beta/MAPK11 protein (**ab117219**)

Lane 2 : Recombinant Human p38 gamma/MAPK12 protein (**ab117221**)

Lane 3 : Recombinant Human p38 delta/MAPK13 protein (**ab113869**)

Lane 4 : Recombinant Human p38 alpha/MAPK14 protein (**ab82188**)

Predicted band size: 41 kDa



Western blot - Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142)

All lanes : Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142) at 1/1000 dilution

Lane 1 : HeLa cell lysate

Lane 2 : Jurkat cell lysate

Lane 3 : Wild-type HEK-293T cell lysate

Lane 4 : MAPK14 knockout HEK-293T cell lysate

Lysates/proteins at 20 µg per lane.


Performed under reducing conditions.

Predicted band size: 41 kDa

Lanes 1 - 4: Merged signal (red and green). Green - ab32142 observed at 40 kDa. Red - loading control, **ab130007** observed at 125 kDa.

ab32142 was shown to react with p38 in wild-type HEK-293T cells.

Loss of signal was observed when knockout cell line **ab255406** (knockout cell lysate **ab263787**) was used. Wild-type and p38 knockout samples were subjected to SDS-PAGE. ab32142 and Anti-Vinculin antibody [VIN-54] (**ab130007**) were incubated overnight at 4°C at 1 in 1000 dilution and 1 in 20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed (**ab216773**) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed (**ab216776**) secondary antibodies at 1 in 20000 dilution for 1 hour at room temperature before imaging.



All lanes : Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142)

Lane 1 : Wild-type HAP1 cell lysate

Lane 2 : p38 knockout HAP1 cell lysate

Lane 3 : HeLa cell lysate


Lane 4 : Jurkat cell lysate

Lysates/proteins at 20 µg per lane.

Predicted band size: 41 kDa

Lanes 1 - 4: Merged signal (red and green). Green - ab32142 observed at 40 kDa. Red - loading control, **ab8245**, observed at 37 kDa.


ab32142 was shown to specifically react with p38 when p38 knockout samples were used. Wild-type and p38 knockout samples were subjected to SDS-PAGE. ab32142 and **ab8245** (loading control to GAPDH) were diluted 1/1000 and 1/2000 respectively and incubated overnight at 4°C. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed (**ab216773**) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed (**ab216776**) secondary antibodies at 1/10 000 dilution for 1 h at room temperature before imaging.



Immunocytochemistry/ Immunofluorescence - Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142)

Immunocytochemistry/Immunofluorescence analysis of MCF-7 (human breast carcinoma) cells labelling p38 (green) with purified ab32142 at 1/250. Cells were fixed with 4% Paraformaldehyde and permeabilized with 0.1% Triton X-100. **ab150077**, Alexa Fluor[®] 488-conjugated goat anti-rabbit IgG (1/1000) was used as the secondary antibody. Nuclei were counterstained with DAPI (blue).


Secondary Only Control: PBS was used instead of the primary antibody as the negative control.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142)

Ab32142, at a 1/100 dilution, staining p38 in paraffin embedded human skin carcinoma tissue by immunohistochemistry.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.



Western blot - Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142)

All lanes : Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142)


- Lane 1 :** Wild-type HAP1 cell lysate
- Lane 2 :** p38 knockout HAP1 cell lysate
- Lane 3 :** HeLa cell lysate
- Lane 4 :** Jurkat cell lysate

Lysates/proteins at 20 µg per lane.

Predicted band size: 41 kDa

Lanes 1 - 4: Merged signal (red and green). Green - ab32142 observed at 40 kDa. Red - loading control, **ab8245**, observed at 37 kDa.

This western blot image is a comparison between ab32142 and a competitor's top cited rabbit polyclonal antibody.



Western blot - Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142)

All lanes : Anti-p38 beta/MAPK11 + p38 alpha/MAPK14 antibody [Y122] (ab32142) at 1/50000 dilution

- Lane 1 :** HeLa cell lysate
- Lane 2 :** Jurkat cell lysate
- Lane 3 :** K562 cell lysate
- Lane 4 :** MCF7 cell lysate

Predicted band size: 41 kDa

Observed band size: 42 kDa

10 ug protein per lane

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-p38 beta/MAPK11 + p38 alpha/MAPK14
antibody [Y122] (ab32142)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors