


Anti-p50 dynamitin antibody ab124285

画像数 2

製品の概要

製品名	Anti-p50 dynamitin antibody
製品の詳細	Rabbit polyclonal to p50 dynamitin
由来種	Rabbit
アプリケーション	適用あり: WB, IP
種交差性	交差種: Mouse, Human 交差が予測される動物種: Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Dog, Xenopus laevis, Chimpanzee, Rhesus monkey, Xenopus tropicalis, Elephant 
免疫原	Synthetic peptide corresponding to Human p50 dynamitin aa 1-50 (N terminal). Database link: NP_006391.1
ポジティブ・コントロール	Jurkat, 293T, HeLa and NIH3T3 whole cell lysates.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7 Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate
精製度	pH 7 to 8 Immunogen affinity purified
特記事項 (精製)	ab124285 was affinity purified using an epitope specific to p50 dynamitin immobilized on solid support.
ポリ/モノ	ポリクローナル

アイソタイプIgG

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab124285の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/2000 - 1/10000.
IP		Use at 2-10 µg/mg of lysate.

ターゲット情報

機能	Modulates cytoplasmic dynein binding to an organelle, and plays a role in prometaphase chromosome alignment and spindle organization during mitosis. Involved in anchoring microtubules to centrosomes. May play a role in synapse formation during brain development.
配列類似性	Belongs to the dynactin subunit 2 family.
細胞内局在	Cytoplasm > cytoskeleton > centrosome. Membrane.

画像




All lanes : Anti-p50 dynamitin antibody (ab124285) at 0.1 µg/ml

- Lane 1 :** Jurkat whole cell lysate at 50 µg
- Lane 2 :** Jurkat whole cell lysate at 15 µg
- Lane 3 :** 293T whole cell lysate at 50 µg
- Lane 4 :** HeLa whole cell lysate at 50 µg
- Lane 5 :** NIH3T3 whole cell lysate at 50 µg

Developed using the ECL technique.

Exposure time: 30 seconds



ab124285 at 1 µg/ml staining p50 dynamitin in Jurkat whole cell lysate immunoprecipitated using ab124285 at 6 µg/mg lysate (1 mg/IP; 20% of IP loaded/lane). Lane 2: control IgG. Predicted band size : 44 kDa. Exposure time: 3 seconds

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors