

Anti-PDK1 antibody [2H3AA11] ab110335

KO 評価済

6 References 画像数 5

製品の概要

製品名	Anti-PDK1 antibody [2H3AA11]
製品の詳細	Mouse monoclonal [2H3AA11] to PDK1
由来種	Mouse
特異性	This antibody does not cross-react with other isoforms of human PDK (i.e. PDK2, 3 or 4) and does not cross-react with rat or mouse samples.
アプリケーション	適用あり: ICC/IF, WB, Flow Cyt
種交差性	交差種: Human
免疫原	Recombinant full length protein. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	HeLa cells; Human heart and liver tissue lysates; and HL-60 cells
特記事項	<p>This monoclonal antibody to PDK1 has been knockout validated in Western blot. The expected band for PDK1 was observed in wild type cells and the band was not seen in knockout cells.</p> <p>This antibody clone is manufactured by Abcam. If you require a custom buffer formulation or conjugation for your experiments, please contact orders@abcam.com.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> <p>Product was previously marketed under the MitoSciences sub-brand.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	<p>pH: 7.5</p> <p>Preservative: 0.02% Sodium azide</p>

	Constituent: HEPES buffered saline
精製度	Proprietary Purification
特記事項 (精製)	Produced in vitro using hybridomas grown in serum-free medium, and then purified by biochemical fractionation. Purity >95% by SDS-PAGE.
ポリ/モノ	モノクローナル
クローン名	2H3AA11
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab110335の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 1 µg/ml.
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 49 kDa.
Flow Cyt		Use a concentration of 1 µg/ml. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

ターゲット情報

機能	Inhibits the mitochondrial pyruvate dehydrogenase complex by phosphorylation of the E1 alpha subunit, thus contributing to the regulation of glucose metabolism.
組織特異性	Expressed predominantly in the heart.
配列類似性	Belongs to the PDK/BCKDK protein kinase family. Contains 1 histidine kinase domain.
細胞内局在	Mitochondrion matrix.

画像



Western blot - Anti-PDK1 antibody [2H3AA11]
(ab110335)

Lane 1: Wild-type HAP1 cell lysate (20 µg)


Lane 2: PDK1 knockout HAP1 cell lysate (20 µg)

Lane 3: Human fetal heart lysate (20 µg)

Lane 4: LnCaP cell lysate (20 µg)

Lanes 1 - 4: Merged signal (red and green). Green - ab110335 observed at 46 kDa. Red - loading control, [ab8245](#), observed at 37 kDa.

ab110335 detected the expected band for PDK1 in wild-type cells and the band was not seen in knockout cells. Additional cross-reactive bands were observed in the wild-type and knockout cells. Wild-type and PDK1 knockout samples were subjected to SDS-PAGE. ab110335 and [ab8245](#) (loading control to GAPDH) were diluted 1 µg/mL and 1/10 000 respectively and incubated overnight at 4°C. Blots were developed with Goat anti-Mouse IgG H&L (IRDye® 800CW) preadsorbed ([ab216772](#)) and Goat Anti-Rabbit IgG H&L (IRDye® 680RD) preadsorbed ([ab216777](#)) secondary antibodies at 1/10 000 dilution for 1 hour at room temperature before imaging.



Western blot - Anti-PDK1 antibody [2H3AA11]
(ab110335)

All lanes : Anti-PDK1 antibody [2H3AA11] (ab110335) at 1 µg/ml

Lane 1 : Human heart lysate at 10 µg

Lane 2 : Human liver lysate at 10 µg


Lane 3 : Recombinant Human PDK1 at 0.01 µg

Lane 4 : Recombinant Human PDK2 at 0.01 µg

Lane 5 : Recombinant Human PDK3 at 0.01 µg

Lane 6 : Recombinant Human PDK4 at 0.01 µg

Predicted band size: 49 kDa



Western blot - Anti-PDK1 antibody [2H3AA11] (ab110335)

All lanes : Anti-PDK1 antibody [2H3AA11] (ab110335) at 1 µg/ml

Lane 1 : HeLa cells incubated for 24 hours with vehicle


Lane 2 : HeLa cells incubated for 24 hours with 100 µM cobalt chloride (CoCl₂)

Lane 3 : HeLa cells incubated for 24 hours with 1 mM deferoxamine (DFO)

Lysates/proteins at 30 µg per lane.


Predicted band size: 49 kDa

Treatment with either CoCl₂ or DFO results in increased expression of PDK1, which is known to be transcriptionally upregulated in response to hypoxia.



Flow Cytometry - Anti-PDK1 antibody [2H3AA11] (ab110335)

HL-60 cells were stained with 1 µg/mL PDK1 antibody (ab110335) (blue) or an equal amount of an isotype control antibody (red) and analyzed by flow cytometry.



Immunocytochemistry/ Immunofluorescence - Anti-PDK1 antibody [2H3AA11] (ab110335)

Immunocytochemistry image of PDK1 stained HeLa cells. The Upper image shows cells which were paraformaldehyde fixed (4%, 20 minutes) and Triton X-100 permeabilized (0.1%, 15 minutes). The cells were incubated with the PDK1 antibody (ab110335) at 2 µg/mL for 2 hours at room temperature or over night at 4°C. The secondary antibody was (green) Alexa Fluor® 488 goat anti-mouse IgG (H+L) at a 1/1000 dilution for 1 hour. 10% Goat serum was used as the blocking agent for all blocking steps. The target protein locates to the mitochondrial matrix.

The lower image shows cells co-stained with an antibody against HSP60 (red), an enzyme also located in the mitochondrial matrix. The composite image shows an identical mitochondrial pattern for both antibodies indicated by merged orange color.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors