

PE Anti-CD52 antibody [YTH34.5] ab33082

画像数 1

製品の概要

製品名	PE Anti-CD52 antibody [YTH34.5]
製品の詳細	PE Rat monoclonal [YTH34.5] to CD52
由来種	Rat
標識	PE. Ex: 488nm, Em: 575nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Tissue, cells or virus corresponding to CD52. Human Blood Lymphocytes (whole cells)
ポジティブ・コントロール	human peripheral blood lymphocytes
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA
精製度	Protein A purified
特記事項 (精製)	Purified IgG prepared from tissue culture supernatant.
ポリ/モノ	モノクローナル
クローン名	YTH34.5
アイソタイプ	IgG2b

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab33082の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		

追加情報

FACS: Use 10ul of neat antibody to label 10^6 cells in 100ul

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報


機能

May play a role in carrying and orienting carbohydrate, as well as having a more specific role.

細胞内局在

Cell membrane.

画像



Staining of human peripheral blood lymphocytes with 10ul NEAT
ab33082 for 10^6 cells per 100ul.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors