

PE Anti-NCAM1 antibody [MEM-188] ab18277

3 References

製品の概要

製品名	PE Anti-NCAM1 antibody [MEM-188]
製品の詳細	PE Mouse monoclonal [MEM-188] to NCAM1
由来種	Mouse
標識	PE. Ex: 488nm, Em: 575nm
特異性	Human NCAM. Reacts with the 180 kDa isoform characteristic for leukocytes. Reactivity with other NCAM isoforms has not been tested.
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Mouse, Human, Monkey
免疫原	Tissue, cells or virus corresponding to NCAM1. KG-1 human acute myelogenous leukaemia cell line
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituents: PBS, BSA</p>
特記事項(精製)	Purified by size exclusion chromatography.
ポリ/モノ	モノクローナル
クローン名	MEM-188
アイソタイプ	IgG2a

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab18277の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration. Designed for analysis of human blood cells using 20ul reagent per 100ul whole blood or 10 ⁶ cells in a suspension. ab91363 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

ターゲット情報

機能	This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.
配列類似性	Contains 2 fibronectin type-III domains. Contains 5 Ig-like C2-type (immunoglobulin-like) domains.
細胞内局在	Secreted and Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors