

Anti-Penicillin antibody [P2B9] ab116641

2 References

製品の概要

| | |
|----------|---|
| 製品名 | Anti-Penicillin antibody [P2B9] |
| 製品の詳細 | Mouse monoclonal [P2B9] to Penicillin |
| 由来種 | Mouse |
| 特異性 | ab116641 shows reactivity with Chemical samples. |
| アプリケーション | 適用あり: IP, ELISA, RIA |
| 種交差性 | 交差種: Species independent |
| 免疫原 | Penicillin |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term. |
| バッファー | Preservative: 0.05% Sodium azide Constituents: 0.02% Potassium chloride, 0.79% Sodium chloride, 0.14% Potassium phosphate, 0.1% BSA |
| 精製度 | Protein A purified |
| ポリ/モノ | モノクローナル |
| クローン名 | P2B9 |
| アイソタイプ | IgG1 |

アプリケーション

Abpromise保証は、次のテスト済みアプリケーションにおけるab116641の使用に適用されます

The Abpromise guarantee

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|-------------------------------------|
| IP | | Use at an assay dependent dilution. |
| ELISA | | Use at an assay dependent dilution. |
| RIA | | Use at an assay dependent dilution. |

ターゲット情報

関連性 Penicillin is a group of Beta-lactam antibiotics used in the treatment of bacterial infections caused by susceptible, usually Gram-positive, organisms. β -lactam antibiotics work by inhibiting the formation of peptidoglycan cross-links in the bacterial cell wall, which results in cytolysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors