

### PerCP Anti-CD3 antibody [UCHT1], prediluted ab106215

7 References [画像数 1](#)

#### 製品の概要

製品名	PerCP Anti-CD3 antibody [UCHT1], prediluted
製品の詳細	PerCP Mouse monoclonal [UCHT1] to CD3, prediluted
由来種	Mouse
標識	PerCP. Ex: 482nm, Em: 675nm
特異性	ab106215 recognizes the CD3 antigen of the TCR/CD3 complex on mature Human T cells. It reacts with the epsilon chain of the CD3 complex. HLDA I; WS Code T 3 HLDA III; WS Code T 126 HLDA III; WS Code T 471 HLDA VI; WS Code T 6T-CD3.1
アプリケーション	<b>適用あり:</b> Flow Cyt
種交差性	<b>交差種:</b> Human, Non human primates
免疫原	Tissue, cells or virus corresponding to Human CD3. Human thymocytes followed by Sezary T cells
ポジティブ・コントロール	Flow Cyt: Human blood cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.4 Preservative: 0.097% Sodium azide Constituents: 0.2% BSA, 2.68% PBS
精製度	Size exclusion
ポリ/モノ	モノクローナル
クローン名	UCHT1
アイソタイプ	IgG1

## アプリケーション

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab106215の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 10µl for 10 <sup>6</sup> cells. <b>ab118658</b> - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

## 画像



Flow Cytometry analysis of Human peripheral blood cells labeling CD3 with Anti-CD3 antibody [UCHT1], prediluted (PerCP) (ab106215).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors