

PerCP Anti-D4DDDK tag (Binds to FLAG® tag sequence) antibody [M2] ab117514

医薬用外毒物

製品の概要

製品名	PerCP Anti-D4DDDK tag (Binds to FLAG® tag sequence) antibody [M2]
製品の詳細	PerCP Mouse monoclonal [M2] to D4DDDK tag (Binds to FLAG® tag sequence)
由来種	Mouse
標識	PerCP. Ex: 482nm, Em: 675nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Species independent
免疫原	Synthetic peptide corresponding to D4DDDK tag (Binds to FLAG® tag sequence).
特記事項	<p>FLAG® is a registered trade mark of Sigma Aldrich Biotechnology LP. It is used here for informational purposes only.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Lyophilized:Store lyophilized protein at -20 C. Reconstitute with 1 ml of ddH2O and store at 2-8°C in the dark and use within 6 months.
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.4 Preservative: 0.01% Sodium azide Constituents: 0.58% Sodium chloride, 0.1% BSA, 1.7% Sucrose, 95% Sodium phosphate
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	M2

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab117514の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		1/500 - 1/5000. ab118658 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

ターゲット情報

関連性 This is a useful tool for the localisation and characterisation of DDDDK tagged proteins (Binds to FLAG® tag sequence).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors